Liquid Handling · Easy Handling!

# Dispensette®

F I R S T C L A S S · B R A N D

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**31** Operating Manual

**59** Mode d'emploi

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	Page
<b>Safety Instructions</b>	<b>32</b>
<b>Functions and Limitations of Use</b>	<b>33</b>
<b>Dispenser Selection Chart</b>	<b>36</b>
<b>Operating Elements</b>	<b>37</b>
<b>First Steps</b>	<b>38</b>
<b>Assembly</b>	<b>38</b>
<b>Priming</b>	<b>40</b>
<b>Dispensing</b>	<b>41</b>
<b>Accessories</b>	<b>42</b>
<b>Error Limits (Nominal Volume, Partial Volume)</b>	<b>45</b>
<b>Checking the Volume (Calibration)</b>	<b>46</b>
<b>Adjustment</b>	<b>47</b>
<b>Cleaning</b>	<b>48</b>
<b>Cleaning/Replacing Valves</b>	<b>50</b>
<b>Autoclaving</b>	<b>51</b>
<b>Ordering Information</b>	<b>52</b>
<b>Accessories · Spare Parts</b>	<b>54</b>
<b>Troubleshooting</b>	<b>57</b>
<b>Repairs and Warranty Information</b>	<b>58</b>
<b>Disposal</b>	<b>58</b>

## Safety Instructions

---

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

### Please read the following carefully!

1. Every user must read and understand this operating manual before operation.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves.
3. Observe all specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid the buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
5. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations. Observe operating exclusions (see page 34)! If in doubt, contact the manufacturer or supplier.
6. Always use the instrument in such a way that neither the user nor any other person is endangered. When dispensing, the discharge tube must always point away from you or any other person. Avoid splashes. Only dispense into suitable vessels.
7. Never press down the piston when the discharge tube closure is attached.
8. Never remove the discharge tube or the SafetyPrime™ recirculation valve while the dispensing cylinder is filled.
9. Reagents can accumulate in the screw cap of the discharge tube. Thus, it should be cleaned regularly.
10. For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over.
11. Never carry the mounted instrument by the cylinder sleeve or the valve block. Breakage or loosening of the cylinder may also lead to personal injury from chemicals (see page 39, Fig. 6).
12. Never use force on the instrument. Use smooth gentle movements to operate the piston upwards and downwards.
13. Use only original manufacturer's accessories and spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
14. Always check the instrument for visual damage before use. If there is a sign of a potential malfunction (e.g., piston difficult to move, sticking valves or leakage), immediately stop dispensing. Consult the 'Troubleshooting' section of this manual (see page 57), and contact the manufacturer if needed.

## Functions and Limitations of Use


---

The bottle-top dispenser Dispensette® III and Dispensette® Organic are designed for dispensing liquids directly from the reservoir bottle. They are available in the models: digital, analog-adjustable and fixed-volume.


The instruments are, according to the requirements of the DIN EN ISO 8655-5, conformity certified and optionally equipped with SafetyPrime™ recirculation valve.

## Dispensette® III (red color code)

Digital · Easy Calibration

Analog-adjustable

Fixed-volume

## Dispensette® Organic (yellow color code)

Digital · Easy Calibration

Analog-adjustable

Fixed-volume

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

### Dispensette® III

Borosilicate glass,  $Al_2O_3$ -ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium, PP (screw cap).

### Dispensette® Organic

Borosilicate glass,  $Al_2O_3$ -ceramic, ETFE, FEP, PFA, PTFE, tantalum, PP (screw cap).

If a higher chemical resistance is required, please use a PTFE screw cap and a ETFE/PTFE bottle adapter (Accessories, pages 54-56).

### Note:

For dispensing HF, we recommend the use of the Dispensette® TA bottle-top dispenser with platinum-iridium valve spring. This is described in a separate operating manual.

# Functions and Limitations of Use

## Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- use temperature from +15 °C to +40 °C (from 59 °F to 104 °F) of instrument and reagent
- vapor pressure up to max. 600 mbar. Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling.
- kinematic viscosity up to 500 mm<sup>2</sup>/s  
(dynamic viscosity [mPas] = kinematic viscosity [mm<sup>2</sup>/s] x density [g/cm<sup>3</sup>])
- Density: Dispensette® III / Dispensette® Organic: up to 2.2 g/cm<sup>3</sup>

## Operating Limitations

Liquids, which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or concentrated alkaline solutions). If the piston becomes difficult to move, the instrument should be cleaned immediately. (page 48).

When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

The Dispensette® is designed for general laboratory applications and complies with the relevant standards, e.g. DIN EN ISO 8655. Compatibility of the instrument for a specific application (e.g., trace material analysis, food sector etc.) must be checked by the user. Approvals for specific applications, e.g. for production and administration of food, pharmaceuticals or cosmetics are not available.

## Operating Exclusions

**Dispensette® III** never use with:

- liquids attacking Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide\*)
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- liquids which are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>)
- hydrochloric acid > 20 % and nitric acid > 30 %
- tetrahydrofuran
- trifluoroacetic acid
- explosive liquids (e.g., carbon disulfide)
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument
- liquids attacking PP (screw cap)\*\*

**Dispensette® Organic** never use with:

- liquids attacking Al<sub>2</sub>O<sub>3</sub>-ceramic, tantalum, ETFE, FEP, PFA and PTFE (e.g., dissolved sodium azide\*)
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- bases and saline solutions
- explosive liquids (e.g., carbon disulfide)
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument
- liquids attacking PP (screw cap)\*\*

\* Dissolved sodium azide permitted up to a concentration of max. 0.1%.

\*\* When stronger chemical resistance is needed, use the PTFE screw caps. (Accessories page 54)

### Storage Conditions

Store the instrument and accessories only in cleaned condition in a cool and dry place.

Storage temperature: from  $-20\text{ }^{\circ}\text{C}$  to  $+50\text{ }^{\circ}\text{C}$   
(from  $-4\text{ }^{\circ}\text{F}$  to  $122\text{ }^{\circ}\text{F}$ ).

### Recommended Application Range

**Dispensette® III** (color code red): Its broad range of application permits bottle dispensing of aggressive reagents, including concentrated acids such as  $\text{H}_3\text{PO}_4$ ,  $\text{H}_2\text{SO}_4$ , bases like NaOH, KOH, saline solutions, as well as many organic solvents.

**Dispensette® Organic** (yellow color code) is ideal for dispensing of organic solvents including chlorinated and fluorinated hydrocarbons (e.g., trichlorotrifluoroethane and dichloromethane), concentrated acids (e.g., HCl and  $\text{HNO}_3$ ), trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides.

### Note:

For guidelines on selecting the right dispenser observe the corresponding Operating Exclusions and the "Dispenser selection chart" on the next page.

For dispensing HF, we recommend the use of the Dispensette® TA bottle-top dispenser with platinum-iridium valve spring. This is described in a separate operating manual.

# Dispenser Selection Chart

Reagent	Disp. III	Disp. Organ.
Acetaldehyde	+	+
Acetic acid (glacial), 100%	+	+
Acetic acid, 96%	+	+
Acetic anhydride		+
Acetone	+	+
Acetonitrile	+	+
Acetophenone		+
Acetyl chloride		+
Acetylacetone	+	+
Acrylic acid	+	+
Acrylonitrile	+	+
Adipic acid	+	
Allyl alcohol	+	+
Aluminium chloride	+	
Amino acids	+	
Ammonia, 20%	+	+
Ammonia, 20 - 30%		+
Ammonium chloride	+	
Ammonium fluoride	+	
Ammonium sulfate	+	
n-Amyl acetate	+	+
Amyl alcohol (Pentanol)	+	+
Amyl chloride (Chloropentane)		+
Aniline	+	+
Barium chloride	+	
Benzaldehyde	+	+
Benzene (Benzol)	+	+
Benzine (Petroleum benzene) bp 70-180 °C		+
Benzyl chloride	+	+
Benzyl alcohol	+	+
Benzylamine	+	+
Benzylchloride	+	+
Boric acid, 10%	+	+
Bromobenzene	+	+
Bromonaphthalene		+
Butanediol	+	+
1-Butanol	+	+
n-Butyl acetate	+	+
Butyl methyl ether	+	+
Butylamine	+	+
Butyric acid	+	+
Calcium carbonate	+	
Calcium chloride	+	
Calcium hydroxide	+	
Calcium hypochlorite	+	
Carbon tetrachloride		+
Chloro naphthalene	+	+
Chloroacetaldehyde, 45%	+	+
Chloroacetic acid	+	+
Chloroacetone	+	+
Chlorobenzene	+	+
Chlorobutane	+	+
Chloroform		+
Chlorosulfonic acid		+
Chromic acid, 50%	+	+
Chromosulfuric acid	+	
Copper sulfate	+	
Cresol		+
Cumene (Isopropyl benzene)	+	+

Reagent	Disp. III	Disp. Organ.
Cyclohexane		+
Cyclohexanone	+	+
Cyclopentane		+
Decane	+	+
1-Decanol	+	+
Dibenzyl ether	+	+
Dichloroacetic acid		+
Dichlorobenzene	+	+
Dichloroethane		+
Dichloroethylene		+
Dichloromethane	+	+
Diesel oil (Heating oil), bp 250-350 °C		+
Diethanolamine	+	+
Diethyl ether		+
Diethylamine	+	+
1,2 Diethylbenzene	+	+
Diethylene glycol	+	+
Dimethyl sulfoxide (DMSO)	+	+
Dimethylaniline	+	
Dimethylformamide (DMF)	+	+
1,4 Dioxane		+
Diphenyl ether	+	+
Essential oil		+
Ethanol	+	+
Ethanolamine	+	+
Ethyl acetate	+	+
Ethylbenzene		+
Ethylene chloride	+	
Fluoroacetic acid		+
Formaldehyde, 40%	+	
Formamide	+	+
Formic acid, 100%		+
Glycerol	+	+
Glycol (Ethylene glycol)	+	+
Glycolic acid, 50%	+	
Heating oil (Diesel oil), bp 250-350 °C		+
Heptane		+
Hexane		+
Hexanoic acid	+	+
Hexanol	+	+
Hydriodic acid	+	+
Hydrobromic acid		+
Hydrochloric acid, 20%	+	+
Hydrochloric acid, 20-37%		+
Hydrogen peroxide, 35%		+
Isoamyl alcohol	+	+
Isobutanol	+	+
Isooctane		+
Isopropanol (2-Propanol)	+	+
Isopropyl ether	+	+
Lactic acid	+	
Methanol	+	+
Methoxybenzene	+	+
Methyl benzoate	+	+
Methyl butyl ether	+	+
Methyl ethyl ketone	+	+
Methyl formate	+	+
Methyl propyl ketone	+	+

Reagent	Disp. III	Disp. Organ.
Methylene chloride		+
Mineral oil (Engine oil)	+	+
Monochloroacetic acid	+	+
Nitric acid, 30%	+	+
Nitric acid, 30-70%*		+
Nitrobenzene	+	+
Oleic acid	+	+
Oxalic acid	+	
n-Pentane		+
Peracetic acid		+
Perchloric acid	+	+
Perchloroethylene		+
Petroleum, bp 180-220 °C		+
Petroleum ether, bp 40-70 °C		+
Phenol	+	+
Phenylethanol	+	+
Phenylhydrazine	+	+
Phosphoric acid, 85%	+	+
Phosphoric acid, 85% + Sulfuric acid, 98%, 1:1		+
Piperidine	+	+
Potassium chloride	+	
Potassium dichromate	+	
Potassium hydroxide	+	
Potassium permanganate	+	
Propionic acid	+	+
Propylene glycol (Propanediol)	+	+
Pyridine	+	+
Pyruvic acid	+	+
Salicylaldehyde	+	+
Scintillation fluid	+	+
Silver acetate	+	
Silver nitrate	+	
Sodium acetate	+	
Sodium chloride	+	
Sodium dichromate	+	
Sodium fluoride	+	
Sodium hydroxide, 30%	+	
Sodium hypochlorite	+	
Sulfuric acid, 98%	+	+
Tartaric acid	+	
Tetrachloroethylene		+
Tetrahydrofuran (THF) **		+
Tetramethylammonium hydroxide	+	
Toluene		+
Trichloroacetic acid		+
Trichlorobenzene		+
Trichloroethane		+
Trichloroethylene		+
Trichlorotrifluoro ethane		+
Triethanolamine	+	+
Triethylene glycol	+	+
Trifluoro ethane	+	+
Trifluoroacetic acid (TFA)	+	+
Turpentine		+
Urea	+	
Xylene		+
Zinc chloride, 10%	+	
Zinc sulfate, 10%	+	


\* use ETFE/PTFE bottle adapter

\*\* use PTFE seal

**Hydrofluoric acid (HF): for dispensing HF, we recommend the use of the Dispensette® TA bottle-top dispenser with platinum-iridium valve spring.**

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BrandTech, Inc. Status as of: 0914/12

# Operating Elements

Filling seal for autoclaving application

Filling and recirculation tube

Mounting tool

## Is everything in the package?

Confirm that your package includes:

Bottle-top dispenser Dispensette®, discharge tube, telescoping filling tube, SafetyPrime™ recirculation valve and recirculation tube (optional), mounting tool, different bottle adapters, O-ring FKM (for autoclaving), a performance certificate and this operating manual.

Nominal volume, ml	Adapters bottle thread, PP	Filling tube Length, mm
0,5	24-25, 28/ S 28, 32-33, 38 mm	125-240
1, 2, 5, 10	28/ S 28, 32-33, 38 mm	125-240
25, 50, 100	32-33, 38 mm	170-330

## Assembly

### Warning:

Wear protective clothing, eye protection and gloves! Follow all safety instructions and observe limitations of use and operating limitations (page 32-34).


### 1. Check sealing washer(s)

Before mounting the SafetyPrime™ recirculation valve or the discharge tube make sure that the sealing washer is inserted.

### 2. Mounting the SafetyPrime™ recirculation valve (optional)

Push the SafetyPrime™ recirculation valve approx. 2 mm into the discharge tube and firmly finger-tighten the locking nut (Fig. 2). Check the SafetyPrime™ recirculation valve for a tight fit.


### 3. Mounting the Discharge tube

Push the discharge tube (optional with SafetyPrime™ recirculation valve, Fig. 3') approx. 2 mm into the valve block and firmly finger-tighten the locking nut (Fig. 3). Check the discharge tube for a tight fit.


### Note:

Recirculation valve and discharge tube must correspond to the model. Pay attention to the color code. On the next day firmly tighten up the locking nut.

## 4. Mounting the filling tube and recirculation tube

Adjust the length of the telescoping filling tube to the bottle height and attach it. If a SafetyPrime™ recirculation valve (optional) is used, the recirculation tube must also be mounted. Insert the recirculation tube with the opening pointing outward (Fig. 4).

## 5. Mounting and aligning the instrument on a bottle

Screw the instrument (GL 45 thread) onto the reagent bottle and then align the discharge tube with the bottle label. By turning the valve block. (Fig. 5). To prevent tipping over use a bottle stand for small bottles.

### Note:

For bottles with other thread sizes, select a suitable adapter.

### Dispensette® III and Dispensette® Organic:

The adapters supplied with the instrument are made of polypropylene (PP), and can only be used for media which do not attack PP.

If a higher chemical resistance is required, please use a EFTE/PTFE bottle adapter (Accessories, page 54).

### Warning:

Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids. When mounted to a reagent bottle, always carry the instrument as shown in the figure (Fig. 6) and always store it in the upright position!

## Warning:

Never press down the piston when the screw cap is screwed on!  
Avoid splashing the reagent! The reagent can drip out from the discharge tube and screw cap.

## Note:

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Avoid splashes.

### Instruments with SafetyPrime™ recirculation valve:

1. Open the screw cap of the dispensing tube (Fig. 1). For safety, hold the discharge tube orifice on the inner wall of a suitable receiving vessel.
2. Set valve to 'Recirculate' (Fig. 2).
3. For priming gently pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure 5 times (Fig. 3).
4. Turn valve to 'Dispense' (Fig. 4).
5. To avoid splashes when priming hold the discharge tube on the inner wall of a suitable receiving vessel and dispense liquid to prime the discharge tube until it is bubble-free. Wipe away any remaining drops from the discharge tube. (Fig. 5).

### Instruments without SafetyPrime™ recirculation valve:

1. Open the screw cap of the discharge tube (see instrument with SafetyPrime™ recirculation valve Fig. 1).  
To avoid splashes, hold discharge tube orifice on the inner wall of a suitable receiving vessel.
2. For priming pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure approximately 5 times until the discharge tube is bubble-free (Fig. 6).

## 1. Setting the volume

**Digital:** Rotate the volume-setting wheel until the desired volume is indicated (mechanical counter).

**Analog-adjustable:** Loosen the volume selector thumb screw one-half turn (1), set the pointer to the desired volume (2) and then retighten the volume thumb screw (3).

**Fixed-volume:** The volume is non-adjustable and cannot be changed.

## 2. Dispensing

### Warning!

Wear protective clothing, eye protection and gloves! Liquid may accumulate in the screw cap. To avoid splashes dispense slowly. Follow all safety instructions and observe limitations of use and operating limitations (page 32-33).

- Remove screw cap from the discharge tube.
- When using instruments equipped with the SafetyPrime™ recirculation valve, turn the valve to "Dispensing".
- Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.
- Gently lift the piston until the upper stop and then depress piston slowly and steadily with minimal force until the lower stop (Fig. d).
- Wipe off the discharge tube against the inner wall of the receiving vessel.
- Reattach screw cap to discharge tube (Fig. f).

### Attention:

After use, always leave the piston in the down position.

## Accessories

The following optional accessories are available for the bottle-top dispenser Dispensette®:

### SafetyPrime™ recirculation valve

The SafetyPrime™ recirculation valve (see accessories, page 55) enables priming of the instrument without loss of medium. Always use the designated recirculation valve that corresponds to the instrument model. For assembly, see "Assembly" on page 10 (Fig. 2).

### Flexible discharge tube

For serial dispensing the flexible discharge tube can be used (Accessories page 55). The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops.

The coil of the tubing can be stretched to a length of the 800 mm max..The entire coil must lie in regular loops and must not be twisted.

The parts in contact with the media are made of:  
Borosilicate glass, Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, PTFE, platinum-iridium.  
Never use for:

- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Peroxides, as they are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>).

Additionally the Operating Exclusions of the instrument apply.

**For mounting,** attach the tube holder onto the valve block (Fig. a) and mount the receiver tube. Slide the dispensing tube with the flexible discharge tube approx. 2 mm further on the valve block, and tighten the lock nut hand-tight. Use a bottle stand (Fig. b) (see accessories, page 56).

### Warning:

There should be no visible damage to the discharge tube (e.g. kinks or the like). Each time you are going to use the tubing, examine it carefully! To dispense aggressive liquids, you should take safety measures in addition to the normal precautions. We recommend use of a protective shield. The bottle must be supported using a bottle stand. To help avoid reagent splashing from the tube, always grip the tube firmly by the handle and replace into the holder after use. For cleaning rinse the tube carefully. Do not dismantle!


### Drying tube

Use of a drying tube, filled with a suitable absorbent (purchased separately), might be necessary for moisture- and CO<sub>2</sub>- sensitive media ('Accessories', page 56).

**For mounting,** unscrew the air vent cap (Fig. a) and screw the filled drying tube in (Fig. b). Place the PTFE sealing ring on the bottle thread (Fig. c) and screw the instrument onto the bottle.

#### Note:

If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

### Sealing ring for valve block

For highly volatile media we recommend sealing the connection from the valve bloc to the bottle with the PTFE sealing ring (see accessories, page 56).

**For mounting,** place the PTFE sealing ring on the bottle thread or the installed adapter (Fig. c) and screw the instrument onto the bottle.

### Air vent cap for micro filter with Luer-cone

For sterile media we recommend the air vent cap with Luer-cone to attach a micro filter. This provides increased protection against contamination by displacement air (see accessories, page 56).

**For mounting,** unscrew the air vent cap (Fig. a) screw in the air vent cap with a Luer cone (Fig. d). Place the PTFE sealing ring on the bottle thread (Fig. c) and screw the instrument onto the bottle. Insert a commercially available sterile filter into the Luer cone (Fig. e).

### Discharge tube with Luer-Lock attachment for micro filter

The dispensing tube with a Luer Lock attachment enables the connection of a microfilter for sterile filtration. The parts which come in contact with the medium are:

Borosilicate glass, Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium and PP (Luer-Lock attachment).

Never use for:

- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- Peroxides, as they are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>).

Additionally observe the Operating Exclusions of the instrument and micro filter.

Make sure that the sealing washer is inserted.

**For mounting**, push the dispensing tube with the Luer Lock (optionally with the SafetyPrime™ recirculation valve) approx. 2 mm into the valve block, and firmly finger-tighten the locking nut (page 38, Figs. 3 and 3'). Check the discharge tube (eventually the SafetyPrime™ recirculation valve) for a tight fit. A commercially available sterile filter can be mounted onto the Luer Lock connector.

#### Note:

Please follow general instructions when handling sterile media. The increased flow resistance can lead to liquid leaking at the upper edge of the dispensing cylinder. To keep any leaking of liquid to a minimum, we recommend using gentle force when dispensing and the use of a filter with a large filter surface.

Error limits related to the nominal capacity (= maximum volume) indicated on the instrument, obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F). Testing takes place according DIN EN ISO 8655-6 with a completely filled instrument and with uniform and smooth dispensing.

Type Digital • Easy Calibration is manufactured under U.S. Patent 5,957,330.

## Error limits Dispensette®

Nominal volume ml	A* ≤ ±		CV* ≤	
	%	µl	%	µl
0,5	1.0	5	0.2	1
1	0.5	5	0.1	1
2	0.5	10	0.1	2
5	0.5	25	0.1	5
10	0.5	50	0.1	10
25	0.5	125	0.1	25
50	0.5	250	0.1	50
100	0.5	500	0.1	100

\* A = Accuracy, CV = Coefficient of Variation

## Partial volume

The percentage values for A and CV are relative to the nominal volume ( $V_N$ ) and must be converted for partial volumes ( $V_T$ ).

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

e.g.	volume	A* ≤ ±		CV* ≤	
		%	µl	%	µl
$V_N$	25.0	0.5	125	0.1	25
$V_T = 50\% N$	12.5	1.0	125	0.2	25
$V_T = 10\% N$	2.5	5.0	125	1.0	25

\* A = Accuracy, CV = Coefficient of Variation

### Note:

The error limits are well within the limits of DIN EN ISO 8655-5. The maximum error limit for a single measurement can be calculated  $EL = A + 2 CV$  (e.g. for volume 25 ml:  $125 \mu l + 2 \times 25 \mu l = 175 \mu l$ ).

# Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to correspond with individual requirements. The complete testing procedure (SOP) can be downloaded at [www.brandtech.com](http://www.brandtech.com). In addition, you can also carry out function checks at shorter intervals, e.g. dispensing the nominal volume into a volumetric test flask (volumetric flask with 3 marks, DAkkS calibrated). For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from [www.brandtech.com](http://www.brandtech.com).

Gravimetric volume testing according to DIN EN ISO 8655-6 (for measurement conditions, see 'Error Limits', page 45) is performed as follows:

## 1. Preparation of the instrument

Clean the instrument ('Cleaning', page 48-51), fill it with distilled H<sub>2</sub>O and then prime it carefully.

## 2. Check the volume

- 10 dispensing operations with distilled H<sub>2</sub>O in 3 Volume ranges (100 %, 50 %, 10 %) are recommended.
- For filling pull up the piston gently until the upper stop of the volume set.
- For discharge depress piston slowly and steadily without force until the lower stop.
- Wipe off the tip of discharge tube.
- Weigh the dispensed quantity on an analytical balance. (Please follow the operating manual of the balance manufacturer.)
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

## 3. Calculations

### Mean volume

$x_i$  = results of weighings  
 $n$  = number of weighings

$Z$  = correction factor  
(e. g., 1.0029 µl/mg at 20 °C, 1013 hPa)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

### Accuracy

$$A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

$V_0$  = nominal volume

### Standard deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

### Coefficient of variation

$$CV\% = \frac{100 \cdot s}{\bar{V}}$$

After a long period of usage an adjustment of the instrument might be necessary.

- Calibrate for example at nominal volume (see page 46).
- Calculate mean volume (result of weighing) (see page 46).
- Adjust the instrument (to the calculated mean volume).
- After the adjustment, further calibration is necessary to confirm appropriate adjustment.


### Example:

Gravimetric testing yields a delivered volume of 9.90 ml with a set volume of 10 ml.

### Type Digital

1. Open housing by sliding the latch to the left and removing the front (Fig. 1).
2. Lift gear lock lever to release. With this action, the adjustment cover plate comes off (Fig. 2).
3. Pull the red knob to disengage the gears and set the display to actual delivered volume (e.g., 9.90 ml) (Fig. 3).
4. Reposition red knob and gear lock lever to their original positions (Fig. 4).
5. Replace housing and slide the latch to the right (Fig. 5). Alteration of factory setting is indicated by a red recalibration flag (Fig. 6).


### Digital

### Type Analog-adjustable

1. Insert the pin of the mounting tool into the cover plate (Fig. 1), and break it off with a rotating motion (Fig. 2).
2. Insert the pin of the mounting tool into the adjustment screw (Fig. 3) and rotate to the left in order to increase the dispensing volume, or rotate to the right to decrease the dispensing volume (e.g. for an actual value of 9.97 ml, rotate approx. 1/2 turn to the left).
3. The change in the adjustment is indicated by an exposed red ring (Fig. 4).

### Type Analog-adjustable

## Adjustment range

Nominal volume	Digital max. +/-	Analog/Fix max. +/-	One rotation corresponds to
0.5 ml	-	5 $\mu$ l	~ 3 $\mu$ l
1 ml	-	6 $\mu$ l	~ 15 $\mu$ l
2 ml	24 $\mu$ l	12 $\mu$ l	~ 15 $\mu$ l
5 ml	60 $\mu$ l	30 $\mu$ l	~ 35 $\mu$ l
10 ml	120 $\mu$ l	60 $\mu$ l	~ 65 $\mu$ l
25 ml	300 $\mu$ l	150 $\mu$ l	~ 130 $\mu$ l
50 ml	600 $\mu$ l	300 $\mu$ l	~ 265 $\mu$ l
100 ml	-	600 $\mu$ l	~ 400 $\mu$ l

# Cleaning

The instrument must be cleaned in the following situations to assure correct operation:

- immediately when the piston is difficult to move
- before changing the reagent
- prior to long term storage
- prior to dismantling the instrument
- prior to autoclaving
- prior to changing the valve
- regularly when using liquids which form deposits (e.g., crystallizing liquids)
- regularly when liquids accumulate in the screw cap


## Warning!

The cylinder, valves, telescoping filling tube and discharge tube contain reagent! Never remove the discharge tube or the SafetyPrime™ recirculation valve while the dispensing cylinder is filled. Point the valves and tube openings away from your body. Wear protective clothing, eye protection and appropriate hand protection.

## Cleaning

For proper cleaning and removal of any deposits, the piston must always be completely pulled out of the cylinder after rinsing.

1. Screw the instrument onto an empty bottle and empty it completely by dispensing (Fig. 1). If the instrument is equipped with SafetyPrime™ recirculation valve, it must be emptied in the "dispensing" and "recirculating" setting.
2. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.
3. If the instrument is equipped with a SafetyPrime™ recirculation valve, after rinsing the instrument, it must also be rinsed in the "recirculating" setting (Fig. 3).
4. Pull out the recirculation tube and the telescoping filling tube.

## Note:

Never exchange pistons between instruments!


## 5. Loosen piston.

### a) Analog-adjusted and fixed-volume model

Hold the housing securely and unscrew the piston completely by turning it to the left.


Carefully pull out the piston.

**Do not remove the housing!**

### b) Digital type

Carry out assembly and dismantling at the maximum volume setting only.



Move the latch to the left and remove the cylinder sleeve.

Place the tip of the mounting tool in the receiving hole, and turn the mounting tool counter-clockwise to loosen the piston mount. Then carefully pull out the piston.

6. Clean piston and cylinder with a bottle-brush (Analog and fixed type see Fig. 6a, Digital type see Fig. 6b). If necessary carefully remove deposits at the edge of the glass cylinder.

7. Then flush all parts of the instrument with deionized water.

8. Insert the piston completely into the cylinder and then reassemble the instrument.

**Note:**

### Digital type

The red stop segment must engage the cylinder's stop ring.


Right!

Wrong!

## Cleaning/replacing the valves

### 1. Filling valve

- Pull out the recirculation tube and the telescoping filling tube (Fig. a).
- Use the mounting tool to unscrew the filling valve (Fig. b).
- If the sealing ring is contaminated or damaged, carefully remove it with a pair of curved forceps (Fig. c).
- Insert cleaned or new sealing ring.
- Screw in the cleaned (e.g., in an ultrasonic bath) or new filling valve first by hand and then tighten it with the mounting tool.

### 2. Discharge valve

The discharge valve is incorporated in the discharge tube. For cleaning see page 48. If necessary dismantle discharge tube and clean it e.g. in an ultrasonic bath. Mount cleaned or new discharge tube see page 38.

### 3. SafetyPrime™ recirculation valve

For cleaning see page 48. If necessary dismantle recirculation valve and clean it e.g. in an ultrasonic bath. Mount cleaned or new recirculation valve see page 38.

#### Note:

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the ball valve is merely stuck.

In this case, loosen the ball valve using light pressure, for example, with a 200 µl plastic pipette tip (see the figure at the side).


This instrument is autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes according to DIN EN 285.

## 1. Prior to the first autoclaving

Remove built-in inelastic PTFE-sealing ring (R) and replace with the supplied elastic sealing ring made of FKM (S).

### Note:



FKM has limited chemical resistance!

R = seal S = elastic O-ring

## 2. Preparation for autoclaving

- Before autoclaving, the instrument must be cleaned carefully (page 48-50). Lift the piston completely out of the cylinder (page 49). For autoclaving with mounted filling tube, the filling valve with olive-shaped nozzle is recommended (page 55).
- Loosen screw connections (Fig. b).
- Pull out the recirculation tube and the telescoping filling tube.
- Check that the filling valve is securely seated (Fig. d). If the rigid PFTE ring has been installed, the filling valve must also be loosened. In the digital model, also check that the piston mount is securely seated (Fig. d').
- Autoclave the instrument with the piston pushed all the way downward. Avoid contact with metallic surfaces. The mounting tool can also be autoclaved.

### Note:

Do not reassemble the instrument until it has cooled down to room temperature (Cooling time approx. 2 hours).

After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them.

It is the user's responsibility to ensure effective autoclaving.

## Ordering Information

### Dispensette® III, Digital · Easy Calibration

Capacity ml	Subdivision ml	without SafetyPrime™ recirculation valve Cat. No.	with SafetyPrime™ recirculation valve Cat. No.
0.2 - 2	0.01	4701 320	4701 321
0.5 - 5	0.02	4701 330	4701 331
1 - 10	0.05	4701 340	4701 341
2.5 - 25	0.1	4701 350	4701 351
5 - 50	0.2	4701 360	4701 361

### Dispensette® III, Analog-adjustable

Capacity ml	Subdivision ml	without SafetyPrime™ recirculation valve Cat. No.	with SafetyPrime™ recirculation valve Cat. No.
0.05 - 0.5	0.01	4701 100	4701 101
0.2 - 2	0.05	4701 120	4701 121
0.5 - 5	0.1	4701 130	4701 131
1 - 10	0.2	4701 140	4701 141
2.5 - 25	0.5	4701 150	4701 151
5 - 50	1.0	4701 160	4701 161
10 - 100	1.0	4701 170	4701 171

### Dispensette® III, Fixed-volume

Capacity ml	without SafetyPrime™ recirculation valve Cat. No.	with SafetyPrime™ recirculation valve Cat. No.
1	4701 210	4701 211
2	4701 220	4701 221
5	4701 230	4701 231
10	4701 240	4701 241

#### Note:

Items supplied see page 38.

## Dispensette® Organic, Digital · Easy Calibration

Capacity ml	Subdivision ml	without SafetyPrime™ recirculation valve Cat. No.	with SafetyPrime™ recirculation valve Cat. No.
0.5 - 5	0.02	4731 330	4731 331
1 - 10	0.05	4731 340	4731 341
2.5 - 25	0.1	4731 350	4731 351
5 - 50	0.2	4731 360	4731 361

## Dispensette® Organic, Analog-adjustable

Capacity ml	Subdivision ml	without SafetyPrime™ recirculation valve Cat. No.	with SafetyPrime™ recirculation valve Cat. No.
0.5 - 5	0.1	4731 130	4731 131
1 - 10	0.2	4731 140	4731 141
2.5 - 25	0.5	4731 150	4731 151
5 - 50	1.0	4731 160	4731 161
10 - 100	1.0	4731 170	4731 171

### Note:

For dispensing HF, we recommend the use of the Dispensette® TA bottle-top dispenser with platinum-iridium valve spring. This is described in a separate operating manual.

## Accessories · Spare Parts

**Bottle adapters** PP or ETFE/PTFE. Adapters of ETFE/PTFE offer higher chemical resistance.

Outer thread mm	for bottle thread/ fits ground joint mm	Material	Pack of	Cat. No.
33	24-25	PP	3	27048 21
33	28/ S* 28	PP	3	27048 22
45	32-33	PP	3	27048 28
45	38	PP	3	27048 27
45	S* 40	PP	3	27048 29
set of addapters	24-25, 28, 32-33, 38	PP	of 1 each	27048 70
33	24-25	ETFE	1	7043 75
33	28/ S* 28	ETFE	1	7043 78
45	32-33	ETFE	1	7043 98
45	38	ETFE	1	7043 99
45	S* 40	PTFE	1	7043 91
33	STj 19/32	PP	3	27048 36
33	STj 24/40	PP	3	27048 38
33	STj 29/42	PP	3	27048 40

\* buttress thread


### Discharge tubes with integrated valve

Pack of 1.

Description	Nominal volume ml	Shape	Length mm	Cat. No.
for Dispensette® III	0.5, 1, 2, 5, 10	fine tip	90	7079 15
	5, 10	standard	90	7079 16
	25, 50, 100	standard	120	7079 17
	25, 50, 100	fine tip	120	7079 18
for Dispensette® Organic	0.5, 1, 2, 5, 10	fine tip	90	7079 35
	5, 10	standard	90	7079 36
	25, 50, 100	standard	120	7079 37
	25, 50, 100	fine tip	120	7079 38

### Telescoping filling tubes

FEP. Adjusts to various bottle heights. Pack of 1.

Nominal volume ml	Outer Ø mm	Length mm	Cat. No.
0.5, 1, 2, 5, 10	6	70-140	7042 02
		125-240	7042 03
		195-350	7042 08
		250-480	7042 01
25, 50, 100	7,6	170-330	7042 04
		250-480	7042 05

### SafetyPrime™ recirculation valves

Pack of 1.

Description	Cat. No.
for Dispensette® III 1-100 ml	7060 80
for Dispensette® III 0.5 ml	7060 81
for Dispensette® Organic	7060 90

### Filling valve with olive-shaped nozzle made of PEEK

For frequent autoclaving with the filling tube mounted, this filling valve with tube nozzle is recommended. Observe limited chemical resistance of PEEK! Pack of 1.

Description	Nominal volume ml	Cat. No.
for Dispensette® III, Dispensette® Organic	0.5, 1, 2, 5, 10	6637
for Dispensette® III, Dispensette® Organic	25, 50, 100	6638

### Flexible discharge tubing

PTFE, coiled, length 800 mm, with safety handle. Pack of 1.

Nominal volume ml	Discharge tube Outer Ø mm	Inner Ø mm	Cat. No.
1, 2, 5, 10	3	2	7079 25*
25, 50, 100	4,5	3	7079 26*

\* not suitable for HF and Peroxide

### Filling valve with sealing washer

Pack of 1.

Description	Nominal volume ml	Cat. No.
for Dispensette® III, Dispensette® Organic	0.5, 1, 2, 5, 10	6697
for Dispensette® III, Dispensette® Organic	25, 50, 100	6698

### Screw cap with fastener

Pack of 1.

Description	Nominal volume ml	Cat. No.
PP, red, for Dispensette® III	0.5, 1, 2, 5, 10	7060 11
	25, 50, 100	7060 12
PP, yellow, for Dispensette® Organic	0.5, 1, 2, 5, 10	7060 13
	25, 50, 100	7060 14
PTFE*, white, for Dispensette® all types	0.5, 1, 2, 5, 10	7060 16
	25, 50, 100	7060 17

\* PTFE, if higher chemical resistance is requested

## Accessories · Spare Parts

### Seals

PTFE. Spare seals for discharge tube, SafetyPrime™ and filling valve.  
Pack of 5 each type.

**Cat. No.** 6696

### Sealing ring for valve block

PTFE, for highly volatile reagents.  
Pack of 1.

**Cat. No.** 7044 86

### Sealing ring for filling valve

FKM (Fluororo elastomer), only for autoclaving applications  
Pack of 5.


**Cat. No.** 6694

### Drying tube

Drying tube and seal, without drying agent.  
Pack of 1.


**Cat. No.** 7079 30

### Air vent cap for micro filter with Luer-cone

PP. Air vent cap and PTFE-sealing ring.  
Pack of 1 each.

**Cat. No.** 7044 95

### Recirculation tube

Pack of 1.

**Cat. No.** 8317

**Discharge tube with Luer-Lock attachment for micro filter**  
FEP/PP. Pack of 1.

**Cat. No.** 7079 28\*

### Calibrating-, mounting-tool

Pack of 1.

**Cat. No.** 6687

\* not suitable for HF and Peroxide

### Bottle Stand

PP. Support rod 325 mm, Base plate 220 x 160 mm.  
Pack of 1

**Cat. No.** 7042 75

Problem	Possible cause	Corrective action
Piston difficult to move	Formation of crystals, dirty	Stop dispensing immediately. Loosen piston with circular motion, but do not disassemble. Follow all cleaning instructions (page 48-50).
Filling not possible	Volume adjusted to minimum setting	Set to required volume (see page 41).
	Filling valve stuck	Clean the filling valve. If the valve is stuck use a 200 µl pipette tip to loosen it (see page 50). If necessary replace the filling valve with sealing washer.
Dispensing not possible	Discharge valve stuck	Clean discharge valve. If necessary replace discharge tube. (see page 50).
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Seal not inserted, knurled locking nuts not firmly connected	Double check that seal is in place, and that the knurled locking nuts at the discharge tube and the SafetyPrime™ recirculation valve are firmly seated and finger tight.
	The instrument has not been primed	Prime the instrument (see page 41).
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged	Cleaning procedure (see page 48-50). Tighten the valves using the mounting tool. If necessary, replace the valves and sealing washers.
	Recirculation tube not connected	Connect recirculation tube (see page 38, Fig. 3).
Dispensed volume is too low	Discharge tube is loose	Push the discharge tube on firmly.
	Filling tube is loose or damaged	Cleaning procedure (see page 48-50). Push the filling tube on firmly. If necessary, cut off approx. 1 cm of the tube at the upper end and re-connect it or replace filling tube (see page 50).
	Filling valve is loose or damaged	Cleaning procedure (see page 48-50). Tighten the valves using the mounting tool. If necessary, replace filling valves and sealing washers.
Leaking liquid between instrument and bottle	Recirculation tube not connected	Connect recirculation tube (see page 38, Fig. 3).
	Volatile reagent dispensed without sealing ring for valve block	Mount sealing ring for valve block (see page 43).

## Repairs, Warranty and Disposal

---

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

**For safety reasons, instruments returned for checks and repairs must be clean and decontaminated!**

### Return for Repair

**Important!** Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

Therefore contact BrandTech Scientific, Inc. and obtain return authorization **before** sending your instrument for service.

Return the instrument, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number. Include an exact description of the type of malfunction and the media used.

## Warranty

---

We shall not be liable for the consequences of improper handling, use, servicing or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used. For length of warranty period please see our warranty card enclosed with the product.

## Disposal

---

For the disposal of instruments, please observe the relevant national disposal regulations.

Subject to technical modification without notice. Errors excepted.

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BrandTech Scientific, Inc. no acepta aparatos que no hayan sido debidamente limpiados y descontaminados.
BrandTech Scientific, Inc. · 11 Bokum Road · Essex, CT 06426- 1506 · USA