
**FINNPIPETTE COLOUR
MULTICHANNEL
FINNPIPETTE MULTICANAUX
DIGITAL-COLORS**
**INSTRUCTIONS FOR USE
GUIDE D'UTILISATION**

Thermo
ELECTRON CORPORATION

Thermo Electron Oy
P.O.Box 100, Fin-01621 Vantaa, Finland
Tel. +358-9-329 100, fax -358-9-3291 0414
www.thermo.com/finnpipette

ENGLISH

PRODUCT DESCRIPTION	3
PIPETTE OPERATION	3
PIPETTING TECHNIQUES	4
MAINTENANCE	5
CALIBRATION	6
PACKAGE	8
TROUBLE SHOOTING	8
SPARE PARTS	9-10
FIGURES	18

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PRODUCT DESCRIPTION

FINNPIPETTE COLOUR MULTICHANNELS are handheld, adjustable volume, general purpose micropipettes for sampling and dispensing of accurate liquid volumes. They operate on the air displacement principle (i.e. an air interface) and use detachable, disposable tips or tip bands.

FINNPIPETTE COLOUR MULTICHANNELS are available as either 8- 12-channel models with a tip ejector system. Each model can be supplied in either of two volume ranges covering a total of 5 to 300 µl. The two volume ranges are:

Order No	Channel	Volume rang	Scale	Increment
4173200	8	5-50 µl	0,5 µl	
4173210	8	50-300 µl	5,0 µl	
4173220	12	5-50 µl	0,5 µl	
4173230	12	50-300 µl	5,0 µl	

PIPETTE MATERIALS

FINNPIPETTE COLOUR MULTICHANNELS are made of mechanically durable and chemically resistant materials.

Operating Cap	PP	Tip ejector	PA
Operating Button	PA	Tip cones	PC
Handle	PA	Pipette body	PA
Tip ejector pusher	PA	Tip ejector plate	PC

DESCRIPTION OF TIPS

FinnTips are recommended for use with the FINNPIPETTE COLOUR MULTICHANNEL.

They are made of natural colour polypropylene, generally regarded as the only contamination free material suitable for tips. FinnTips are also autoclavable (121°C).

Following accessories are useful with handling of pipettes and pipetting of reagents:

Stand	9420290
Shelf hanger	2206040

DIGITAL DISPLAY (Fig. 1)

The adjusted delivery volume is clearly indicated by a digital display on the pipette handle.

PIPETTE OPERATION

SETTING THE DELIVERY VOLUME (Fig. 2)

The delivery volume is set using the operating button on the top of the pipette. To increase the delivery volume turn the operating button counterclockwise. To decrease turn it clockwise.

Ensure that the desired delivery volume clicks into place and that the digits for the volume are completely visible on the handle display.

Do not use force to try to set volumes outside of the pipette's specified volume range!

TIP EJECTION (Fig. 3)

Each pipette is fitted with a tip ejector system to help eliminate the possibility of contamination. For tip ejection, two methods are available:

- A. If one-handed operation is preferred, twist the pipette handle so that the tip ejector pusher and tip ejector on the body assembly are in line. Press the tip ejector pusher with your thumb directing the tips towards a suitable waste receptable.
- B. In two-handed operation, press the tip ejector on the body assembly directing the tips towards a suitable waste receptacle. With some tip types this may be an easier and lighter method for user if the tips are very tight.

PIPETTING TECHNIQUES**Figures 4-6:**

- A = Ready position
B = First stop
C = Second stop

FORWARD TECHNIQUE (Fig. 4)

Ensure that clean tips are firmly pushed on to the tip cones of the pipette.

The operation of the pipette is controlled by the operating button. For the best possible accuracy:

Operate the button slowly and with steady speed at all times, particularly when working with liquids having high viscosity. Never let the operating button snap back. When working with serum or biological fluids prime the tips with the liquid before pipetting. This is done by filling and emptying the tips.

Hold the pipette vertically (10 degrees maximum from vertical) during the liquid intake. For maximum hand control, hold the pipette tightly in the hand with the "Grippy" resting on the index finger.

Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button to the first stop.
2. Dip the tips just under the surface of the liquid in the reservoir (2-3 mm) and slowly release the operating button. Withdraw the tips from the liquid touching them against the edge the reservoir to remove excess liquid.
3. Deliver the liquid by gently depressing the operating button to the first stop. After a delay of about one second, continue to depress the operating button all the way down to the second stop. This action will empty the tips.
4. Release the operating button to the ready position.

If necessary, change the tips and continue with the pipetting.

REVERSE TECHNIQUE (Fig. 5)

The reverse technique is suitable for dispensing liquids having high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button all the way down to the second stop.
2. Dip the tips just under the surface of the liquid in the reservoir (2-3 mm) and slowly release the operating button.
This action will fill the tips. Withdraw the tips from the liquid touching them against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the operating button to the first stop. Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.
4. The remaining liquid is either discarded with the tips or pipetted back into the container.

REPETITIVE TECHNIQUE (Fig. 6)

The repetitive technique offers a rapid and simple procedure for repeated delivery of the same volume. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button all the way down to the second stop.
2. Dip the tips just under the surface of the liquid in the reservoir (2-3 mm) and slowly release the operating button.
This action will fill the tips. Withdraw the tips from the liquid touching them against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the operating button to the first stop. Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.
4. Continue pipetting by repeating the steps 2 and 3.

MAINTENANCE

When the FINNPIPETTE COLOUR MULTICHANNEL is not in use make sure it is safely stored in a vertical position, preferably on a Finnpipette Stand.

SHORT TERM CHECKING

The pipette should be checked at the beginning of each day for dust and dirt on outside surfaces of the pipette. Particular attention should be paid to the tip cones. Solvents should not be used for cleaning the pipette.

LONG TERM MAINTENANCE

If the pipette is used daily it should be checked at least twice a year. The servicing procedure is as follows:

1. Unscrew the two screws on the tip ejector plate and remove it.
2. Unscrew the eight screws on the body assembly cover and remove it.
3. Remove the tip ejector with the spring.
4. Depress the operating button all the way down. Remove the tip cone holder with the tip cones.

5. Clean the pistons and O-rings with a dry napless cloth and grease them slightly with grease (included in the package).

6. Reassemble the pipette.

In the event of a malfunction clean the pipette and send it to the Thermo Electron representative.

Please, enclose a complete list of liquids which have entered the pipette.

CALIBRATION

All Finnpipettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. The pipettes are constructed to permit re-adjustment for liquids of different temperature and viscosity.

DEVICE REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range	readable graduation
under 10 µl	0.001 mg
10-100 µl	0.01 mg
above 100 µl	0.1 mg

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ($\pm 0.5^{\circ}\text{C}$) temperature of water, pipette and air between 15°C to 30°C .

The relative humidity must be above 50%. Especially with volumes under 50 µl the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

Procedure to check calibration

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume. A new tip is first pre-wetted 3-5 times and a series of ten pipettings is done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume.

Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the inaccuracy (A) and imprecision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the calculated results are within the selected limits, the adjustment of the pipette is correct.

TABLE 1: Maximum permissible errors according ISO8655

Range	Channel	Volume µl	Inaccuracy µl	%	Imprecision s.d. µl	Imprecision cv%
5-50 µl	8	50	±0.75	±1.5	0.35	0.7
		5	±0.25	±5.0	0.10	2.0
50-300 µl	8	300	±3.0	±1.0	0.9	0.3
		50	±2.3	±4.6	0.8	1.5
5-50 µl	12	50	±0.75	±1.5	0.35	0.7
		5	±0.25	±5.0	0.10	2.0
50-300 µl	12	300	±3.0	±1.0	0.9	0.3
		50	±2.3	±4.6	0.8	1.5

ADJUSTMENT (Fig. 7)

Adjustment is done with the service tool.

1. Place the service tool into the openings of the calibration nut at the top of the handle.
2. Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.
3. After adjustment check the calibration according to the instructions above.

FORMULAS FOR CALCULATING RESULTS

Conversion of mass to volume

$$V = (w + e) \times Z$$

V = volume (µl)
 w = weight (mg)
 e = evaporation loss (mg)
 Z = conversion factor for mg/µl conversion

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6 mg = 0.2 mg/s).

Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2 mg (10 s x 0.2 mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary.

The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032 µl/mg at 22°C and 95 kPa. See the conversion table on page 18.

Inaccuracy (systematic error)

Inaccuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0$$

A = inaccuracy
 \bar{V} = mean volume
 V_0 = nominal volume

Inaccuracy can be expressed as a relative value:

$$A\% = 100\% \times A / V_0$$

Imprecision (random error)

Imprecision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$S = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

s = standards deviation
 \bar{V} = mean volume
 n = number of measurements

Standard deviation can be expressed as a relative value (CV)

$$CV = 100\% \times S / \bar{V}$$

PACKAGE

The FINNPIPETTE COLOUR MULTICHANNELS are shipped in a custom made, cardboard package.

The package contains the following items:

1. FINNPIPETTEMCP
2. Instruction Manual
3. Tube of Grease
4. Service tool
5. Samples of Tips

TROUBLE SHOOTING


The following table is a guide to possible problems and their solutions.

Fault	Possible reason	Solution
Leakage	Tip incorrectly attached	Attach firmly
	Foreign bodies between tips and tip cones	Clean tip cones thoroughly and attach new tips
	Foreign bodies between O-rings and pistons	Clean and grease O-rings and the pistons in the cylinders. Use grease
	Insufficient amount of grease on pistons and O-rings	Grease accordingly
	O-rings damaged	Change the O-rings
Inaccurate dispensing	Incorrect operation	Follow instructions carefully
	Tips incorrectly attached	Attach firmly
	Calibration altered: e.g. caused by misuse	Recalibrate according to instructions
Inaccurate dispensing with certain liquids	Unsuitable calibration. Many high viscosity liquids require recalibration	Recalibrate with the liquids in question

CAUTION!

The Finnpipette is designed to allow easy in-lab service. If you would prefer to have us or your local representative service your pipette, please make sure that the pipette has been decontaminated before you send it to us. Please note that the postal authorities in your country may prohibit or restrict the shipment of contaminated material by mail.

Description			Code
1. Handle (incl. item 18)		5-50 µl 50-300 µl	2205130 2205140
2. Tip ejector		5-50 µl 50-300 µl	10590140 10590150
3. Tip ejector spring			1130111
4. Body assembly cover screw			0202061
5. Tip ejector screw			0202021
6. Tip ejector plate	8-channel 8-channel 12-channel 12-channel	5-50 µl 50-300 µl 5-50 µl 50-300 µl	10590300 10590310 10590320 10590330
7. Body assembly cover	8-channel 8-channel 12-channel 12-channel	5-50 µl 50-300 µl 5-50 µl 50-300 µl	10590240 10590250 10590260 10590270
8. Tip cone		5-50 µl 50-300 µl	10590340 10590350
9. O-ring		5-50 µl 50-300 µl	1030480 1030260
10. Collar	8-12-channel 8-12-channel	5-50 µl 50-300 µl	1056030 1056040
11. Tip cone holder	8-channel 12-channel 8-channel 12-channel	5-50 µl 5-50 µl 50-300 µl 50-300 µl	1056090 1056100 1056120 1056130
12. Label			1527110
13. Tip ejector pusher holder			1054820
14. Spring support			1056630
15. Tip ejector pusher spring			1130570
16. Tip ejector pusher		5-50 µl 50-300 µl	10588776 10588786
17. Service tool			10589070
18. Cap			10588850
19. Body	8-channel 8-channel 12-channel 12-channel	5-50 µl 50-300 µl 5-50 µl 50-300 µl	2205090 2205100 2205110 2205120
20. Piston holder	8-channel 12-channel		10590770 10590810
21. Holder spring	8-channel 12-channel		1130680 1130720
22. Locking pieces			1058180
23. Piston		5-50 µl 50-300 µl	2201240 2202940
24. Piston spring	8-12-channel		1131630

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Finnpipette® and Finntip® are registered trademarks of Thermo Electron Oy.