Gilson Pipetting Aid

FRANÇAIS	4 – 11
ENGLISH	12 – 19
DEUTSCH	20 – 27
ESPAÑOL	28 – 35
РУССКИЙ	36 – 43

PRINTED IN POLAND

LT801478B/E/JULY 2003 PRINTED IN POLAND


World Headquarters
Gilson, Inc.
3000 W. Beltline Hwy., P.O. Box 620027, Middleton, WI 53562-0027, USA
Telephone: (1) 800-445-7661 or (1) 608-836-1551 • Fax: (1) 608-831-4451

Gilson SAS
19 avenue des Entrepreneurs, 95400 Villiers-le-Bel, France
Telephone: (33) 1-34-29-50-00 • Fax: (33) 1-34-29-50-20

Gilson Pipetting Aid

1

FRANÇAIS

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TABLE OF CONTENTS

1 - INTRODUCTION	13
2 - WORK SAFETY INSTRUCTIONS	13
3 - LIMITATIONS OF USE	14
4 - PREPARING FOR USE	14
5 - ASPIRATING AND DISPENSING LIQUIDS	14
6 - TROUBLESHOOTING	15
7 - REPLACING THE FILTER AND CLEANING THE VALVES	17
8 - CHARGING THE BATTERY	17
9 - MAINTENANCE	18
10 - COMPONENTS	18
11 - ORDERING INFORMATION	19
12 - REPLACEMENT PARTS	19

1 - INTRODUCTION

The **Gilson Pipetting Aid** is a device intended for pipetting liquids with the use of measuring pipettes. It can work with all types of glass or plastic pipettes in the volume range from 1 mL to 100 mL.

The **Gilson Pipetting Aid** is equipped with replaceable membrane filters, which protect the mechanism of the device from contamination by the vapours of solution being drawn into the pipette.

The two-step drawing speed control system enables both very fast dispensing of large volumes and precise measuring of small volumes.

Fig. 1 shows the external parts of the pipetting aid with a description of the materials used.

2 - WORK SAFETY INSTRUCTIONS

Before starting to work with the **Gilson Pipetting Aid** every user should read these operating instructions carefully.

Using the device inconsistently with the operating instructions may result in damaging the device.

The device should be maintained only by your Gilson Distributor, otherwise the manufacturer will be relieved from any liability under the warranty.

- When working with the **Gilson Pipetting Aid** general safety regulations regarding risks related to laboratory work should be observed. Protective clothing, goggles and gloves should be worn.
- **Gilson Pipetting Aid** shall be used only for measuring liquids in conditions specified by the manufacturer, which are limited due to the chemical and mechanical resistance of the device, as well as the user safety.
- **Gilson Pipetting Aid** should not be used where the risk of an explosion is present.
- The information and instructions provided by the manufacturers of the reagents must be observed.
- Only original spare parts and accessories (See Chapter 12 for user replacement parts), recommended by the manufacturer, shall be used.
- Only the original charger, supplied by the manufacturer, shall be used for charging the battery.
- In case of incorrect functioning of the pipetting aid, work shall be stopped. The device shall be cleaned in accordance with the operating instructions and/or sent for repair to your **Gilson** Distributor.
- In the case of mechanical damage to the casing, the device shall be immediately sent for repair to an authorised service centre.
- The use of excessive force during work shall be avoided.
- If you need to dispose of the battery, for any reason, please do it in an environmentally friendly way, in accordance with the applicable regulations.

3 - LIMITATIONS OF USE

- The **Gilson Pipetting Aid** may not be used for measuring substances, the vapours of which damage the following plastics: PP, SI, ABS, EPDM, POM.
- Gilson Pipetting Aid** may not be used where the risk of an explosion risk is present.
- Flammable liquids shall not be measured - in particular substances with flash- point below 0 °C (ether, acetone).
- Gilson Pipetting Aid** shall not be used for drawing acids with a concentration above 1 mol.
- Gilson Pipetting Aid** shall not be used for drawing solutions with a temperature above 50 °C.
- Gilson Pipetting Aid** may work in temperature range from +10 °C to +35 °C.

4 - PREPARING FOR USE

Before starting, check that the battery has not discharged. To check this, press the aspiration button (fig. 1A) and watch to see if the low battery indicator lights up (fig. 1E). If the indicator does not light up, the pipetting aid is ready for use.

If the low battery indicator lights up, the pipetting aid shall be left for charging (recommended action), or the charger shall be connected and the work continued during charging.

The **Gilson Pipetting Aid** may be charged only with the original charger.

The mains voltage shall conform with the specification on the charger.

Charging shall be done in accordance with Chapter 8 of this instruction manual.

5 - ASPIRATING AND DISPENSING LIQUIDS

Attaching a pipette

Before attaching a pipette, check that the pipette is not damaged, has no dents or sharp edges in the part to be gripped. Check that the part to be gripped is dry.

The pipette shall be gripped as close to the upper end as possible and carefully inserted into the pipette holder until resistance is noticed (fig. 2A). Do not apply excessive force so as not to damage thin pipettes and to avoid injury risk. A pipette that has been correctly attached and sealed in the holder should not tilt to the side.

After attaching a pipette, hold the pipetting aid so as to keep the pipette in the vertical position. The shape of the casing makes it possible to lay down the pipetting aid with a pipette attached. After use, it is recommended not to leave the pipetting aid with a pipette attached for a long period, for example overnight or over a weekend.

Do not lay the pipetting aid down, when there is liquid in the pipette.

Filling the pipette

Before aspirating, set the suction speed using the HIGH/LOW switch (fig. 1C):

- HIGH position - fast aspirating
- LOW position - slow aspirating

The LOW position is recommended for pipettes with a volume up to 5 mL, and the HIGH position for pipettes with a volume greater than 5 mL. Hold the pipetting aid so that the pipette is in the vertical position, immerse the pipette end in the liquid to be aspirated (fig. 2B), and press the aspiration button gently. The pipette filling speed depends on how far the aspiration button is pressed. The further the button is pressed, the faster the liquid is aspirated into the pipette.

Aspirating a slightly greater volume than required is recommended (due to the meniscus above the required volume mark). Adjust the aspiration speed, particularly in the final filling stage, so as not to overfill the pipette.

Setting the volume

After filling, dry the end of the pipette with absorbent paper (one that does not leave impurities), so as to remove any liquid from the outside surface of the pipette. Then set the required liquid volume precisely. Press the dispense button gently (fig. 2C), dispense the excess liquid from the pipette until the meniscus of the liquid aligns exactly with the mark on the pipette that indicates the volume to be dispensed.

Dispensing

Hold the recipient vessel in inclined position, place the pipette end in contact with the inside wall, and press the dispense button gently (fig. 2C). The dispensing intensity may be adjusted depending on how far the dispense button is pressed. The deeper the button is pressed the faster the outflow of liquid from the pipette.

The **Gilson Pipetting Aid** has two dispense modes. The dispense mode is selected with the GRAV/ BLOW switch (fig. 1D).

- GRAV position - the dispense is effected in gravity mode, which means that the liquid flows out of the pipette by its own weight. The outflow speed is adjusted depending on how far the dispense button is pressed.
- BLOW position - the dispense is effected as in the GRAV mode, but when the dispense button is pressed to the end position, the pump is started and the dispense is completed by blowing the liquid from the pipette.

Warning!

In the GRAV mode, the pipette is not completely emptied because of the characteristics of the pipettes used with the pipetting aid.

6 - TROUBLESHOOTING

If the pipetting aid does not function correctly, you should determine the cause and remove the fault, before continuing to use the pipetting aid.

Problem	Possible cause	Action
The pipette falls out (the holding force of the pipette is too small), or tilts to the side too much.	The pipette holder is dirty or wet (fig.1G).	Take the pipette out of the holder. Disassemble the pipette holder and check for mechanical damage. Clean, wash and dry the pipette holder.
	The pipette holder is damaged.	If the holder is damaged, replace it with a new one.
The pump is working, but the pipetting aid does not draw liquid or draws liquid very slowly.	The filter is dirty (fig. 1H).	Take the pipette out of the holder. Disassemble the pipette holder. Take out the filter - check it for impurities - if it is dirty, replace it with a new one.
	The check valve is blocked (fig. 1J).	Check the valve for flooding. If it is wet - rinse it and dry. Check the valve for blockage or displacement due to excessive force during pipette installation - if it is blocked or displaced try to set it in the right position.
Discharged batteries.		Connect the charger to charge the batteries.
	The pipette holder and/or the connector gasket is damaged (fig. 1K).	Check the pipette holder and the connector gasket for mechanical damage - if there is a damage, replace the damaged part with a new one.
Liquid leaks from the pipette (the aspiration and the dispense buttons are not pressed)	The pipette is damaged.	Check the pipette for damage - cracks, dents - if present, replace the pipette with a new one.
	The pipette is inserted incorrectly.	Check that the pipette has been correctly inserted in the holder.
	The pipette holder, the filter or the connector gasket is installed incorrectly.	Check that all parts are present and correctly installed.

The pipette holder and/or the connector gasket is damaged (fig.1G, fig 1K). Check the pipette holder and the connector gasket for mechanical damage - if there is damage, replace the damaged part with a new one.

If the actions specified in the table do not solve the problem, pipetting aid should be sent to your **Gilson** Distributor.

Warning! Before sending the product should be cleaned and decontaminated.

7 - REPLACING THE FILTER AND CLEANING THE VALVES

Note:

The work safety instructions given in Chapter 2 shall be observed when disassembling the pipetting aid.

Replace the filter, if the aspiration efficiency deteriorates. This is because the filter may be dirty after a long period of use.

Note: Only use the original manufacturers filters.

The procedure is shown in (fig. 3).

- Remove the pipette.
- Unscrew the nosepiece (fig. 3A).
- Remove the membrane filter (fig. 3A) and the pipette holder (fig. 3B).
- Rinse the holder with the safety valve using a wash bottle (fig. 3C).
- Blow liquid out of the holder and leave it until it is completely dry.
- Install new membrane filter (fig. 3D) and assemble the device in reverse order.

8 - CHARGING THE BATTERY

The **Gilson Pipetting Aid** has a low battery indicator. If the control lamp (fig. 1E) lights up, when the aspiration button is pressed, it means that the battery is discharged and should be charged.

Note:

The Gilson Pipetting Aid may be charged only with the original charger. The mains voltage must conform with the specification on the charger.

Using other chargers than the original one may damage the battery of the pipetting aid.

Battery charging is regulated by a time circuit, which controls the entire process. When the battery has been charged, the charging circuit disconnects automatically. This is indicated by the control lamp (fig. 1E), which shines continuously.

Charging:

1. Charging temperature: 10 °C to 35 °C.
2. Charging is done with the charger included with each pipetting aid. Battery charging is signalled by the control lamp, which flashes during charging.
3. Full charging time: 11- 14 hours.

4. The charging condition is signalled by the control lamp, which shines continuously when the charging is complete.

The pipetting aid has a NiMH type battery with a capacity of 1700 mAh.

The service life of the battery is about 1000 charging cycles, if used correctly. It is not possible to overcharge the battery, provided that all of the manufacturers instructions are followed.

Note:

Charging should not be interrupted. If the charging process is interrupted, the next charging process should take place after the battery has been completely discharged.

9 - MAINTENANCE

Cleaning and Decontamination

The **Gilson Pipetting Aid** does not require any maintenance. The external parts of the pipetting aid may be cleaned with a swab moistened with isopropanol.

The nosepiece and the pipette holder with the check valve may be autoclaved at 121 °C for 20 minutes. The filters included in the set may be sterilised by autoclaving at 121 °C for not more than 15 minutes.

After autoclaving and cooling the pipette holder, bend it and squeeze it several times in the area of the safety valve.

Storage

The **Gilson Pipetting Aid** shall be stored in a dry place. The allowable storage temperature: - 20 °C to + 50 °C.

When not in use, the pipetting aid shall be placed on the wall nosepiece hanger supplied with the product.

The nosepiece hanger may be fitted on a wall on a hook or stuck with the adhesive tape attached to the hanger. The surface for the hanger should be smooth, clean and degreased. Before sticking the nosepiece hanger, remove the protective film from the tapes attached to the hanger.

Note:

Do not store the pipetting aid with a filled pipette.

10 - COMPONENTS

The pack includes the following components:

1 Gilson Pipetting Aid

1 Battery charger

1 Nosepiece hanger

2 Membrane filters (0.2 µm)

1 Instruction manual

11 - ORDERING INFORMATION

Version	Voltage	Reference number
Gilson Pipetting Aid with EU power supply	220 V / 50 Hz	F161250

* If you want to order a Gilson Pipetting Aid with a US, UK or Australia power supply, please contact your Gilson distributor.

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