Liquid Handling · Easy Handling!

HandyStep® S

F I R S T C L A S S · B R A N D

HandyStep® S

4 - 5

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21 Operating Manual

39 Mode d'emploi

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| | Page |
|------------------------------------------------|-----------|
| Safety Instructions | 23 |
| Function and Limitations of Use | 24 |
| Operating and Control Elements | 25 |
| Inserting the PD-Tip | 26 |
| Setting the Volume | 27 |
| Filling the PD-Tip | 28 |
| Repetitive Dispensing | 29 |
| Ejecting the PD-Tip | 30 |
| Checking the Volume | 31 |
| Accuracy Table | 33 |
| Servicing · Cleaning · UV sterilization | 34 |
| Ordering Information · Accessories | 35 |
| Troubleshooting | 36 |
| Repairs | 37 |
| Warranty Information | 38 |
| Disposal | 38 |

Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves.
When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe the reagent manufacturers' information.
4. Only use the instrument for dispensing liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 24).
If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Do not attempt to make any technical alterations. Do not dismantle the instrument.
9. Before use, check the instrument for visible damages. If there is a sign of a potential malfunction, immediately stop dispensing. Consult the 'Troubleshooting' section of this manual (see page 36), and contact the manufacturer if needed.

Function and Limitations of Use

The HandyStep® S is a repetitive pipette for quick and simple repetitive dispensing of liquids. Volumes from 2 µl to 5 ml can be dispensed with the highest precision and accuracy through combinations of the PD-Tip size used and adjustment of the stroke setting on the stroke-setting wheel.

Limitations of Use

This instrument is designed for dispensing samples, observing the following limits:

- +15 °C to +40 °C (59 °F to 104 °F)
(of instrument and reagents – other temperatures on request)
- vapor pressure up to 500 mbar
- viscosity: 20 mPa s with 50 ml PD-Tips
260 mPa s with 5 ml PD-Tips
977 mPa s with 1 ml PD-Tips

Operating Exclusions

During proper operation, the dispensed sample only comes in contact with the tip and not with the HandyStep® S.

The user has to ensure the compatibility of the instrument with the intended application. In case of third-party dispenser tips, the user must check their compatibility and perform a function and volume test before use.

The instrument should not be used for dispensing liquids that attack polypropylene (e.g., PD-Tip cylinder), polyethylene (PD-Tip piston), LCP (PD-Tip piston 0.1 ml), or PC/PBT and polycarbonate.

Note:

Instrument and tips are not autoclavable (sterile PD-Tips, see page 35). PD-Tips are disposables. For disinfection of the instrument you may use commercially available disinfecting solutions as long as they do not attack PC/PBT and PC.

Operating and Control Elements

| Setting | Top size (ml) | | | |
|---------|---------------|-----|-----|-----|
| | 0.1 | 0.5 | 1 | 2.5 |
| 1 | 3 | 10 | 20 | 50 |
| 2 | 3 | 15 | 30 | 75 |
| 3 | 4 | 20 | 40 | 100 |
| 4 | 5 | 25 | 50 | 125 |
| 5 | 6 | 30 | 60 | 150 |
| 6 | 7 | 35 | 70 | 175 |
| 7 | 8 | 40 | 80 | 200 |
| 8 | 9 | 45 | 90 | 225 |
| 9 | 10 | 50 | 100 | 250 |

Volume (µl)

Backside: volume table

Inserting the PD-Tip

- Use of the 25 ml and 50 ml PD-Tips requires the reusable adapter supplied with the PD-Tips. This is coupled to the PD-Tips with a bayonet lock, and can be detached after use.
- PD-Tips are disposables.

1. Push the locking/
filling lever
to the lower
stop.

2. Swing out the
locking/filling
lever.

3. Insert PD-Tip
straight and
upright from
below.

4. Swing back
the locking/
filling lever.

Adapter

Setting the Volume

The possible dispensing volumes of the HandyStep® S are listed in tables. These are found beneath the transparent clip on the backside of the instrument.

There are 9 half-integral stroke setting positions, from 1 to 5.

1. Look for the desired dispensing volume (Volume) in the table.
2. Use the stroke setting wheel to set to the proper stroke number (Setting) for the volume that corresponds to the PD-Tip (Tip size).
3. Dispensing steps: depending on the PD-Tip capacity (ml), the same desired dispensing volume per step can be delivered a different number of times (see example below). The smaller the number of steps, the higher the accuracy of the dispensed volume.

Example:

Desired dispensing volume per step: **200 µl**

PD-Tip sizes to use:

2.5 ml = stroke setting **4** = max. **11** dispensing steps

5 ml = stroke setting **2** = max. **24** dispensing steps

10 ml = stroke setting **1** = max. **49** dispensing steps

| Setting | Tip size (ml) | | | | Steps |
|-------------|---------------|-----|-----|-----|-------|
| | 0.1 | 0.5 | 1 | 2.5 | |
| 1 | 2 | 10 | 20 | 50 | 49 |
| . | 3 | 15 | 30 | 75 | 32 |
| 2 | 4 | 20 | 40 | 100 | 24 |
| . | 5 | 25 | 50 | 125 | 19 |
| 3 | 6 | 30 | 60 | 150 | 15 |
| . | 7 | 35 | 70 | 175 | 13 |
| 4 | 8 | 40 | 80 | 200 | 11 |
| . | 9 | 45 | 90 | 225 | 10 |
| 5 | 10 | 50 | 100 | 250 | 9 |
| Volume (µl) | | | | | |

| Setting | Tip size (ml) | | | | Steps |
|-------------|---------------|------|------|------|-------|
| | 5 | 10 | 25 | 50 | |
| 1 | 100 | 200 | 500 | 1000 | 49 |
| . | 150 | 300 | 750 | 1500 | 32 |
| 2 | 200 | 400 | 1000 | 2000 | 24 |
| . | 250 | 500 | 1250 | 2500 | 19 |
| 3 | 300 | 600 | 1500 | 3000 | 15 |
| . | 350 | 700 | 1750 | 3500 | 13 |
| 4 | 400 | 800 | 2000 | 4000 | 11 |
| . | 450 | 900 | 2250 | 4500 | 10 |
| 5 | 500 | 1000 | 2500 | 5000 | 9 |
| Volume (µl) | | | | | |

Important!

The volumes listed in the table are the only options available. Since the first dispensing step must be discarded, the number given in the table is always 1 step less.

Filling the PD-Tip

Aspirating Liquid

1. Push the locking/filling lever to the lower stop.
2. Immerse the orifice of the PD-Tip vertically into the liquid to a depth of 3 - 10 mm.
3. Raise the locking/filling lever slowly to the upper stop to avoid the formation of air bubbles.

Note:

Any small air bubble near the piston can be ignored since a locking mechanism prevents the residual volume from being dispensed inadvertently after the final complete dispensing step.


Dispensing Operation

1. Check the volume setting once again.
(Compare with the settings in the volume table.)
2. Wipe off any liquid adhering to the outside of the PD-Tip orifice with a lint-free tissue.
3. Place the orifice of the PD-Tip against the wall of the vessel.
4. Dispense the liquid by pressing down the dispensing lever completely, then let it slide all the way back up.
5. Take care to dispense smoothly and evenly.

Important!

The first dispensing step must be discarded!

Ejecting the PD-Tip

Operation

1. Hold the HandyStep® S over a vessel.
2. Empty the PD-Tip by pushing the locking/filling lever to its lower stop.
3. Swing out the locking/filling lever (the piston will then be decoupled).
4. Hold the HandyStep® S over the waste container. Press the ejection key down, releasing the cylinder. The PD-Tip will be ejected.


Warning!

Residues of the medium might be present on the PD-Tip orifice. The PD-Tip should only be ejected in a manner that does not pose a hazard to the user or other persons.

Note:

The tightness of the seal of the PD-Tips is not guaranteed for repeated dispensing of high viscosity media.

Checking the Volume

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to individual requirements. Gravimetric volume testing according to DIN EN ISO 8655-5 is performed as follows:

1. Preparation of the instrument

Fitting tip. This test can be carried out with PD-Tips of any size. However, size 5 ml is most commonly used for this purpose.

2. Carry out the test

- a) Adjust the stroke setting of HandyStep® S to Step 5 (10 % of nominal volume).
- b) Filling the PD-Tip, immersing it vertically into the testing liquid.
- c) Discard the first step; it only serves to align the internal mechanism.
- d) Place the weighing vessel (containing a small amount of deionized water) on a balance and tare the balance.
- e) Dispense the second step into the weighing vessel. For this, push down the dispensing lever at a steady rate to the stop and hold it there. Then wipe the tip on the vessel wall over a distance of approx. 10 mm.
- f) Enter weighed value into test record.
- g) Repeat points D to F 10 times.
(For the 10th volume test, the 5 ml PD-Tip with a stroke setting of 5 has to be filled a second time – steps A through C.)
- h) Repeat the same testing procedure at stroke settings 3 (6 % of nominal volume) and 1 (2 % of nominal volume).
- i) This results in a total of the 30 weighed values needed.

Checking the Volume

Calculation (for nominal volume)

x_i = Weighing results

n = Number of weighings

Z = Correction factor

(e.g., 1.0029 µl/mg
at 20 °C, 1013 hPa)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

Accuracy*

$$A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

V_0 = Nominal volume

Coefficient of Variation*

$$CV\% = \frac{100}{\bar{V}}$$

Standard Deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation of accuracy (A%) and variation coefficient (CV%):
A% and CV% are calculated according to the formulas for
statistical control.

Note:

Testing instructions (SOPs) and a demo version of the EASYCAL™ 4.0 calibration software are available for download at www.brand.de.

Technical Data

HandyStep® S with PD-Tips from BRAND, 20 °C 'Ex', ^B

| Nominal volume PD-Tip | Volume range (µl) | | A* ≤ ± % | | | CV* ≤ % | | |
|--------------------------|-------------------|------|-----------------------------------------|--------|--------|---------|-----------------------------------------|--------|
| | | | Stroke setting ≈ % of nominal volume | 1 ≈ 2% | 3 ≈ 6% | 5 ≈ 10% | Stroke setting ≈ % of nominal volume | 1 ≈ 2% |
| 0.1 ml | 2 - | 10 | 8.0 | 2.7 | 1.6 | 5.0 | 3.0 | 2.0 |
| 0.5 ml | 10 - | 50 | 4.0 | 1.33 | 0.8 | 1.4 | 0.73 | 0.6 |
| 1.0 ml | 20 - | 100 | 4.0 | 1.33 | 0.8 | 1.0 | 0.38 | 0.4 |
| 1.25 ml | 25 - | 125 | 4.0 | 1.33 | 0.8 | 0.8 | 0.38 | 0.3 |
| 2.5 ml | 50 - | 250 | 3.5 | 1.17 | 0.7 | 0.8 | 0.3 | 0.2 |
| 5.0 ml | 100 - | 500 | 2.5 | 0.83 | 0.5 | 0.6 | 0.27 | 0.2 |
| 10.0 ml | 200 - | 1000 | 1.5 | 0.5 | 0.4 | 0.5 | 0.23 | 0.2 |
| 12.5 ml | 250 - | 1250 | 1.5 | 0.5 | 0.3 | 0.3 | 0.23 | 0.2 |
| 25.0 ml | 500 - | 2500 | 1.5 | 0.5 | 0.3 | 0.4 | 0.23 | 0.2 |
| 50.0 ml | 1000 - | 5000 | 1.5 | 0.5 | 0.3 | 0.4 | 0.23 | 0.15 |

A* = Accuracy, CV* = Coefficient of variation

Error limits refer to the partial volume set relative to the PD-Tip, obtained when instrument, PD-Tip and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. The testing is according to DIN EN ISO 8655-5.

Important!

The HandyStep® S can also be operated with compatible tips from other manufacturers. The user must check their compatibility and perform a function and volume test before use!

Servicing · Cleaning · UV sterilization

HandyStep® S is factory calibrated and maintenance-free.
The instrument must not be dismantled!

In cases of external contamination, we recommend cleaning
the instrument using a wipe with water or isopropanol.

Important!

Avoid liquid getting inside the instrument. In case this happens,
let the instrument be opened and cleaned by an authorized
service center, or send it back to BrandTech.

UV sterilization

The unit can withstand the usual
output of a UV sterilization lamp.
The effects of the UV may cause
some color change.

The rack mount is provided with
an adhesive strip for mounting. It
can be used without the strip in the
Transferpette® S 6-position bench-
top rack.

Ordering Information · Accessories

English

HandyStep® S

pack of 1, 1 shelf/rack mount,
3 PD-Tips (0.1 ml, 1 ml, 10 ml).

Shelf/rack mount
pack of 1.

Cat. No.

27051 10

Cat. No.

7051 30

PD-Tips (Positive Displacement Tips)

| Capacity | Pack of | non sterile Cat. No. | sterile** Cat. No. | BIO-CERT® Cat. No. |
|----------|----------|-------------------------|-----------------------|-----------------------|
| 0.1 ml | 100 | 7024 02 | 7024 04 | 7026 83 |
| 0.5 ml | 100 | 7023 70 | 7023 84 | 7026 84 |
| 1.0 ml | 100 | 7024 06 | 7024 36 | 7026 85 |
| 1.25 ml | 100 | 7023 72 | 7023 86 | 7026 86 |
| 2.5 ml | 100 | 7023 74 | 7023 88 | 7026 88 |
| 5.0 ml | 100 | 7023 76 | 7023 90 | 7026 90 |
| 10.0 ml | 100 | 7024 07 | 7024 38 | 7026 91 |
| 12.5 ml | 100 | 7023 78 | 7023 92 | 7026 92 |
| 25.0 ml* | 50/25*** | 7023 80 | 7023 94 | 7026 94 |
| 50.0 ml* | 25 | 7023 82 | 7023 96 | 7026 96 |

* incl. 1 adapter, ** sterile/free of endotoxin – single wrapped,

*** PD-Tips 25 ml: non sterile pack of 50 / sterile pack of 25

Adapter for 25 and 50 ml

PD-Tips, PP, autoclavable

non sterile, pack of 10

Cat. No.

7023 98

BIO-CERT®, pack of 5

Cat. No.

7026 99

PD-Tip Set

20 PD-Tips each of 0.5, 1.0, 1.25,
2.5, 5.0, 10 and 12.5 ml size

Cat. No.

7023 68

PD-Tips BIO-CERT® are sterile, free of endotoxins, DNA, RNase and ATP.

Troubleshooting

| Problem | Possible cause | Corrective action |
|--------------------------------------------------------------------|---------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| PD-Tip piston cannot be locked after the PD-Tip has been inserted. | Locking/filling lever is not pushed down completely and is not swung out. | First press the ejection key and remove the PD-Tip, then push the locking/filling lever all the way down and swing it out completely. |
| Undefined dispensing volume | The stroke setting wheel is not properly engaged. | Ensure that the stroke setting wheel is securely locked into place. |
| Air bubble beneath the PD-Tip piston becomes larger. | Leaking PD-Tip | Replace PD-Tip |
| A and/or CV exceed tolerances. | Leaking PD-Tip The dispenser tips are difficult to move. | Replace the PD-Tip or use other more pliant dispenser tips. |
| | The dispensing lever has not been pressed completely and evenly. | Press the dispensing lever evenly to stop position. |
| | The first dispensing step was not discarded. | Discard the first dispensing step. |
| The PD-Tip piston cannot be decoupled after emptying completely. | The PD-Tip has been inserted improperly. | Press the ejection key, push down the locking/filling lever completely, and swing it out. |

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

Please note: for the safety of couriers and our employees, and to avoid violation of federal and local laws, only clean instruments free of any chemical, biological or radioactive hazards can be inspected and repaired!

Return for Repair

Important! Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

Therefore contact BrandTech Scientific, Inc. and obtain return authorization **before** sending your instrument for service.

Return the instrument, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number. Include an exact description of the type of malfunction and the media used.

Warranty ---

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used. For length of warranty period please see our warranty card enclosed with the product.

Disposal ---

For the disposal of instruments and PD-Tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice. Errors excepted.

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BrandTech Scientific, Inc. no acepta aparatos que no hayan sido debidamente limpiados y descontaminados.
BrandTech Scientific, Inc. · 11 Bokum Road · Essex, CT 06426-1506 · USA