# PIPETMAN® 8X200

FRANÇAIS


1 - 11

ENGLISH

12 - 22

DEUTSCH

23 - 34


**GILSON®**  
SOLUTIONS AT WORK FOR YOU


# PIPETMAN 8X200




1

2

4

# PIPETMAN 8X200

## SOMMAIRE

---

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### **1. INTRODUCTION**

---

Pipetman 8X200 is an eight-channel Pipetman for liquid measurement and transfer applications requiring the use of 96-well plates. The 8X200 allows continuous volume adjustment from 20 to 200 µl. A digital volumeter indicates the volume setting in microliters, simplifies the volume setting process and reduces the chance of error in reading volumes.

Each nozzle contains its own piston and seal assembly ([Fig. 3](#)). The piston seals are lubricated with silicone grease for low friction and long life. Precision-ground stainless steel pistons contribute to seal integrity and ensure accuracy and precision for each nozzle. Nozzle assemblies are removable and can be individually replaced.

A piston bar, controlled by the operator's thumb, moves the eight pistons simultaneously during operation. Tips are easily ejected by a cam-actuated tip ejector, thus avoiding the need for hand contact with contaminated tips ([Fig. 4](#)).

# PIPETMAN 8X200

The eight-channel liquid end can be rotated to any angle for operator comfort.

The Pipetman 8X200 is a high quality instrument which offers excellent accuracy and precision. The figures given in the following table are only guaranteed when using Gilson tips.

Volume µl	Accuracy (mean error)		Precision (repeatability)	
	Absolute µl	Relative %	Absolute µl	Relative %
20	± 0.5	± 2.5	≤ 0.25	≤ 1.25
50	± 0.5	± 1.0	≤ 0.25	≤ 0.5
100	± 1.0	± 1.0	≤ 0.5	≤ 0.5
200	± 2.0	± 1.0	≤ 1.0	≤ 0.5

These specifications are obtained in forward mode, using a gravimetric method with the temperature of the distilled water and all other conditions stabilized at  $21.5^{\circ} \pm 1.5^{\circ}\text{C}$ . The values given include all components of error due to both normal hand warming and the changing of the tip.

**Performance tests:** Each instrument is inspected and validated according to the Gilson Quality Assurance System. Based on extensive historical data, manufacturing conditions and expertise, and in compliance with ISO standards relative to statistical process control, the assurance level of this instrument performing to specifications is 99.8%. Specifications rely on the quality and consistency for the whole pipetting system; they are guaranteed with Gilson manufactured tips.

## 2. SETTING THE VOLUME

The volumeter consists of three number dials and is read from top to bottom, the bottom figure being the least significant. The three digits indicate the volume selected in microliters. The volume indicated is the volume for each channel.

You can set intermediate volumes using the vernier marks between the digits on the bottom dial. Pipetman 8X200 operates in the range from 20 to 200 µl; the smallest incremental setting is 0.2 µl.

Various volume settings for Pipetman 8X200 are shown opposite:

0	1
6	0
7	8

67  $\mu$ l 108  $\mu$ l

The volume is set by turning either the adjustment screw ([Fig. 1F](#)) or the adjustment knob ([Fig. 2A<sub>2</sub>](#)). When wearing gloves, you may find it easier to use the adjustment knob to set the volume. To obtain maximum accuracy when changing the volumeter setting, follow these recommendations:

- When **decreasing** the volume setting, turn the adjustment ring slowly to reach the required setting, making sure you don't overshoot the mark.
- When **increasing** the volume setting, turn the adjustment ring until you are 1/3 of a turn above the required setting. Then turn the adjustment ring slowly to decrease the volume setting to the required value, making sure you don't overshoot the mark.

### 3. OPERATION

Mount new disposable tips on all eight pipette nozzles ([Fig. 1I](#)) by pressing down on the Gilson Tipack™ rack. Use a back-and-forth or side-to-side rocking motion to ensure that the tips are securely mounted. Press firmly to ensure a positive airtight seal. Set the desired volume using the adjusting screw.

#### ASPIRATION

- 1) Press the push-button ([Fig. 2A](#)) to the first stop ([Fig. 1A](#)). This part of the stroke is the calibrated volume displayed on the digital volume scale.
- 2) Holding the 8X200 vertically, immerse the tips into the samples to an immersion depth of 2-4 mm.
- 3) Allow the push-button to return slowly and smoothly to the top position to aspirate the sample ([Fig. 1B](#)).

# PIPETMAN 8X200

- 4) Wait one second to ensure that the full volume of sample is drawn into each tip.
- 5) Withdraw the tips from the vials. If some liquid remains on the outside of the tips, wipe them carefully with a medical wipe, taking care not to touch the tip orifices.

## DISPENSING

- 1) Dispense the sample by touching the tip ends against the side walls of the receiving wells and pressing the pushbutton slowly to the first stop ([Fig. 1C](#)).
- 2) Wait one second (longer for viscous solutions). Then press the push-button to the second stop (bottom of the stroke, [Fig. 1D](#)) to expell any residual liquid; this is the purge stroke.
- 3) Keeping the push-button fully depressed, withdraw the 8X200 carefully, wiping the tips along the walls of the wells.
- 4) Allow the push-button to return to the uppermost position ([Fig. 1E](#)).
- 5) Discard the tips by twisting the winged tip ejector knob. New tips should be used for each sample to prevent sample carry-over.

## PRE-RINSING

When you fit a new tip (or change the volume to be aspirated), you should pre-rinse the tip. To pre-rinse a tip you aspirate the desired volume and then dispense it back into the same vessel.

The action of pre-rinsing a new tip improves the accuracy and repeatability for the volume of liquid subsequently transferred. This is because when you aspirate a liquid using a pipette a film forms on the inside of the tip. The nature of the film, which is a potential source of error, depends on the liquid being transferred. However, this film remains relatively constant after successive operations with the same tip. Therefore, you pre-rinse the tip to maximize the performance of your pipette.

## PIPETTING LIQUIDS OF VARIOUS DENSITIES

It is possible to compensate for solutions having densities much different to that of water by setting the digital micrometer slightly higher or lower than the required volume. The amount of compensation must be obtained empirically.

For example, when using a 8X200 Pipetman to pipette eight channels of 100 µl of a dense solution, you determine gravimetrically that the volume delivered is actually 85 µl. Try setting the digital micrometer to 115 µl and repeat the measurements. If the volumes delivered are still not sufficiently close to 100 µl, make another slight micrometer adjustment, corresponding to the direction of the error.

 Very dense or viscous liquids, or volatile liquids may not be suitable for air displacement pipetting. Use Gilson Microman® positive displacement pipettes for these liquids.

---

## 4. GILSON DIAMOND TIPS

---

Gilson tips are made to the highest specifications, strict quality control is maintained throughout the manufacturing process:

- They are process clean, meeting today's requirements for dye-free materials (avoiding the possibility of contamination).
- Optimized shape (revised collar for optimum sealing, thin walls, and fine point), making them easier to mount, more flexible, with no vortexing, and improved precision.
- They are free from even microscopic defects, especially at the orifice. All surfaces are smooth and hydrophobic, thereby avoiding the excessive retention of liquids that causes poor accuracy and a lack of precision.
- Mold and cavity references are marked on the collar, ensuring the traceability for quality assurance purpose.
- They form an air-tight seal with the tip-holder, preventing the leaks that cause poor accuracy and a lack of precision.

# PIPETMAN 8X200

 The quality of Gilson Diamond Tips guarantees the accuracy and precision figures specified in **Chapter 1**. The use of inferior quality tips will seriously degrade the performance of your Pipetman.

## D200 Gilson Diamond Tips

These tips are used for volumes between 2 µl and 200 µl and to calibrate Pipetman 8X200.

- They have level marks at 20 and 100 µl which allow you to quickly check the volume.
- They have the Gilson emblem engraved on their collar, ensuring that you have genuine Gilson product.
- They may be autoclavable at 121°C for 20 minutes at 0.1MPa.

---

## 5. RECOMMENDATIONS

---

### PIPETTING GUIDELINES AND PRECAUTIONS

Consistency in all aspects of pipetting procedure will significantly contribute to optimum reproducibility. Pay particular attention to:

- Consistent aspirate/dispense rhythm while pipetting.
- Consistent speed and smoothness of action when pressing and releasing the push-button.
- Consistent push-button pressure at the first stop.
- Consistent immersion depth of all tips.
- Pipetting angle of 20° maximum from the vertical.
- When aspirating, if air bubbles are seen in any of the tips, dispense the liquid into the original vessels, check the tip immersion depth (2-4 mm) and pipette more slowly. If a bubble appears again, eject the tips and use a new set.

## ROTATING THE LIQUID END

It is possible to rotate the eight-channel liquid end for comfort of operation. To do this, first slacken the connecting nut ([Fig. 1G](#))  $\frac{1}{4}$  turn. Rotate the liquid end to the position you want. Re-tighten the connecting nut.

## ACIDS AND CORROSIVE LIQUIDS

Pipetting strong acids or corrosive solutions with Pipetman 8X200 is not recommended. Extensive contact with corrosive fumes may corrode the pistons, result in premature seal wear, and may require premature replacement of the nozzle assemblies ([Fig. 2D](#)).

## TIP EJECTOR

Do not take the eject mechanism apart ([Fig. 1I](#)). If disassembled, reassembly will be impossible.

## AFTER USE

Pipetman 8X200 should be placed on the Gilson Carrousel (reference F144301) when not in use. Do not place the instrument on its side with sample liquid in the tips.

**Pipetman 8X200 is not autoclavable, nor should it be exposed to UV light.**

---

## 6. PERFORMANCE CHECKING

---

### LEAK CONTROL

- 1) Fit eight D200 tips from the Tipack. Use a back-and-forth rocking motion to ensure that the tips are securely attached
- 2) Set the volumeter to 200 and aspirate 200  $\mu\text{l}$  of water.
- 3) Hold the 8X200 vertical and observe the tip ends for 15 seconds. If you observe a drop of water at the end of the tip, or a leak, go to section 7. If no leak is observed, go to the second stage of this procedure.

# PIPETMAN 8X200

## GRAVIMETRIC TEST

This test is based on the determination of the weight of water samples delivered by the pipette. For this test a balance with a sensitivity of  $10^{-5}$  g is needed.

The recommended operating temperature is  $21.5^{\circ} \pm 1.5^{\circ}\text{C}$ , deionized water stabilized at the ambient temperature should be used.

To control evaporation, a small amount of water is placed in the covered weighing vessel.

- 1) Record ambient and water temperatures.
- 2) Fit a tip to channel 1 (extreme left) and set the desired volume by adjusting the volumeter.
- 3) Pre-rinse the tip once. This measurement will not be taken into account for the calculations.
- 4) Cover the receptacle and place it on the balance pan. Reset the electronic balance, or record the value (for a mechanical balance). While waiting for the stabilization, aspirate the sample using the forward mode.
- 5) Open the balance door, retrieve the weighing receptacle, deliver the sample, place the receptacle on the balance pan and close the balance door.
- 6) Allow the display to stabilize, record the value.
- 7) Repeat steps 1 to 6 as many times as required by the test (see [Chapter 6](#), Test Frequency).
- 8) Record the ambient and water temperatures again. Check that the values are still within the recommended limits.
- 9) Convert the results from mg to  $\mu\text{l}$  using the correction factor  $Z^*$ . Use the average of the first and second values as the test temperature to determine the correction factor Z. (For an average temperature of  $21.5^{\circ}\text{C}$  and atmospheric pressure around 1013 hPa, Z factor = 1.0032). Calculate the inaccuracy and imprecision and compare them with the specifications given in the introduction of the manual.

10) Repeat the whole procedure for channel 8 (extreme right).

 \* For more information about Z factor, please refer to the local Gilson representative.

If a problem occurs, or if your Pipetman needs to be recalibrated, send it to your Gilson representative.

## TEST FREQUENCY

- **Functional test:**

Leakage tests should be performed daily.

- **Quick check:**

Perform the “gravimetric test” for four samples, this is usually sufficient to control the accuracy and roughly estimate the reproducibility. It should be performed at least monthly.

- **Routine control:**

To control accuracy and reproducibility, perform the gravimetric test for ten samples. This test should be carried out quarterly.

---

## 7. MAINTENANCE

---

Pipetman 8X200 is solidly constructed to provide exceptional performance and long-term trouble-free service, provided proper pipetting techniques are used.

The piston seal ([Fig. 3C](#)) is designed to withstand extended wear; however, damage to a piston ([Fig. 3B](#)) may impair the seal and cause leakage. If either situation occurs, the damaged nozzle assembly must be replaced.

## TROUBLESHOOTING

Possible causes of leaks or incorrect sampling:

- 1) **Cracked nozzle:** Remove the liquid end cover ([Fig. 1H](#)) and inspect the nozzles ([Fig. 3A](#)). Replace if necessary.

# PIPETMAN 8X200

- 2) **Insufficient lubrication:** Remove the liquid end cover. Lift out the affected nozzle, remove and lubricate the piston with a thin film of silicone grease; then reassemble.
- 3) **Worn piston seal:** Remove the liquid end cover and replace the affected nozzle assemblies.
- 4) **Improper reassembly:** Remove the liquid end cover and check the position of the nozzle assemblies, ensuring that the nozzle locating keys are pressed into the locating slots.
- 5) **Tips do not fit properly:** Use only Gilson D200 Tips and Gilson Tipack racks.

## NOZZLE REMOVAL

Set the pipette volume to zero.

Undo the two Phillips screws situated at the back of the liquid end cover.

Undo the four Phillips screws from the front of the liquid end cover and lift off the front cover (see [Fig. 1H](#)).

Lift out the damaged nozzle assembly ([Fig. 3](#)). This is best done by grasping the nozzle body and pulling up and out from the liquid end. The nozzle body may slide off the piston. If this happens, remove the piston separately.

## NOZZLE REPLACEMENT

In order to secure a nozzle assembly in a station:

- 1) Position the replacement nozzle assembly directly above the vacant position.
- 2) Push the piston clip ([Fig. 3D](#)) into the piston bar slot ([Fig. 2C](#)). The clip will compress slightly as you press it into place.
- 3) Engage the nozzle body locating key ([Fig. 3E](#)) in the locating slot in the liquid end box ([Fig. 2E](#)). This should be a snug fit. The nozzle should not move when you try to twist it. All nozzles should line up evenly.

Replace the front cover and screw the four cover screws back into place. Do not over-tighten these screws or the corners of the cover may crack. Then tighten the two screws in the back of the liquid end cover.

If the problem continues after carrying out nozzle replacement, send the pipette to your Gilson representative.

 **Before returning a pipette, please ensure that it is completely free of any chemical, biological, or radioactive contamination.**

## **8. SPARE PARTS LIST**

You can order the following spare parts through your Gilson representative:

Description	Part Number
Push-button ( <a href="#">Fig. 2A</a> ), set of 5	F1861081
Nozzle/Piston Assembly ( <a href="#">Fig. 3</a> )	F161067
Silicone Grease	F1485012
Spring (Internal)	F1861072
Cover plus screws ( <a href="#">Fig. 1H</a> )	F1861073

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### **Gilson, Inc.**
### **Gilson SAS**