# Transferpette® -8/-12

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**Read carefully before you start!**

**Safety Instructions**

1. Observe general hazard prevention instructions and safety regulations, e. g. wear protective clothing, goggles and gloves.
2. Observe the reagent manufacturer's information.
3. Read this Operating Manual carefully.
4. Always use the instrument in such a way that neither the user nor any other persons are endangered.
5. Only use manufacturer's original parts.

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**Purpose**

The multichannel pipette is an air displacement system for pipetting aqueous solutions with medium density and viscosity into titration plate.

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**Operating Limitations**

The user has to ensure the compatibility of the instrument with the intended application.

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**This instrument cannot be used for:**

- liquids incompatible with polypropylene
- liquids of a very high vapour pressure.

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
Viscous and wetting liquids may impair the accuracy of the volume. The same applies to liquids of a temperature deviating from the ambient temperature by more than ± 5 °C (permissible operating range: 15 - 40 °C).

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
## Handling

### Use proper pipette tips

Pipette performance can only be guaranteed if high-quality pipette tips are used. We consequently recommend the use of PLASTIBRAND® pipette tips.

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## Adjusting the manifold

The **manifold** can be rotated counterclockwise to the most convenient pipetting position.

 counterclockwise rotation for ease of handling

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

**Note:** The manifold must be completely mounted to the handle for accurate pipetting. During the counterclockwise rotation, a clicking noise must be heard.

## Pipetting

Securely attach tips to nose cone assembly. Put your thumb **across** the pipetting key – unlike with conventional pipettes.

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

## Pipetting

### Aspirating the sample

- Rinse the tip once with the sample liquid.
- ① Press the pipetting lever on the side to the **first stop (A)**.
- ② Dip the pipette tip 2-3 millimeters into the sample.

**Note:** Leave the tip in the liquid for about one second longer to avoid the intake of air.

- ③ Let the pipetting lever slide back **slowly**.
- ④ Lightly wipe the tip against the wall of the vessel.

**Note:** Don't lay the instrument horizontal when the tips are full, or liquid may penetrate inside the instrument.

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## Dispensing the sample

Place the pipette tip against the wall of the vessel.

- Press the pipetting key down to the stop (A) and hold it down.

**Note:** For serums and liquids of high-viscosity or low surface tension, wait about 3 seconds to improve the accuracy.

- The blow-out stroke empties the tips completely: Press down to the stop (B).

**Note:** When working at low volumes to improve the accuracy: Rinse the tips with the reagent contained in the vessel.

Wipe the pipette tips against the wall of the vessel.

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


### Ejecting the tips


Press the ejector cap down firmly (see arrow). Tips that are jammed on tightly should be pulled off separately by hand.

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

### After use:

Always store the instrument upright on the **stand** supplied.

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**Note:** To reduce strain on the sealing rings, always store the instrument without the tips attached.

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
### Maintenance


#### Check for tightness

After drawing in the sample: Hold down the instrument vertically for about 10 seconds. If a drop forms at one of the tips,  
➤ "Trouble Shooting", page 19.

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### Lubrication:

After sterilising the manifold several times (more than 5 times, page 17), it is necessary to lubricate the pistons to maintain smooth operation. Dismantle the instrument (➤ "Cleaning", page 15). To lubricate the pistons, use silicone fluid without additives, viscosity of min. 5000 mm<sup>2</sup>/s. A suitable silicone fluid is available as accessory (➤ "Ordering Data", page 27).

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


## Adjustment

The instrument is permanently calibrated for watery solutions. It can be adjusted for solutions of other density and viscosity.

1. Check the volume, determine actual value (page 13).
2. Release lateral closures by turning into position (⊕)
3. Remove ejector cap (A).

4. Use a paper clip to extract the red security plate (B) from the top.
5. Pull out volume setting mechanism (C) to the stop in the direction of its axis (approx. 2 mm).

6. Set to the previously measured actual value. Push back volume setting mechanism (C) to the stop (approx. 2 mm), without turning it.
7. Replace the red security plate.
8. Mount ejector cap (A).
9. Secure ejector cap by turning the lateral closures into position (⊖).

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$$\text{V\%} = \frac{100 s}{\bar{V}}$$

## Checking the volume

**Note:** Detailed Testing Instructions are available under Cat.No. 9976 23.

The gravimetric volume test of the pipette requires the following steps:

(Testing procedures are described e.g. in ISO DIS 8655 Part 3, and DIN 12650 Part 6).

- Adjust volume to nominal capacity.

### Conditioning the pipette:

- To condition the pipette before testing, mount a new pipette tip and pipette five times with testing liquid (H<sub>2</sub>O dist.).  
Discard the pipette tip.

### Carrying out the test (observe instructions by the balance manufacturer):

- Mount a new pipette tip and pre-rinse once with testing liquid.
- Take up testing liquid and pipette into the weighing vessel.
- Weigh the pipetted quantity with an analytical balance.
- Calculate the volume, taking the temperature into account.
- Carry out at least 5 pipetting and weighing operations.

Calculation of accuracy (A%) and coefficient of variation (CV%):  
A and CV are calculated according to the equations used in statistical quality control.

### Calculations (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$x_i$  = results of the weighings  
 $n$  = number of weighings

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
## Components

- ① Handle
- ② Manifold
- ③ Manifold housing
- ④ Nose cone assembly
- ⑤ Spring
- ⑥ Piston
- ⑦ Piston mount
- ⑧ Manifold housing cover

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## Cleaning

### Dismantling the instrument

**Be careful:** Work over a tray to collect any part that may fall out during dismantling.

- Eject the tips.
- Set the instrument to maximum volume.
- Remove the manifold (2) from the handle (1).
- Turn both closures of the manifold housing cover (8) with a coin (rotation through 90°) and slide off the housing.
- Remove both screws which hold the nose cone assembly (4) to the metal support rods of the manifold housing cover (8).
- Remove the nose cone assembly (4), springs (5), and piston mount (7). Do not dismantle the nose cone assembly any further.

### Cleaning the components

- Clean nose cone assembly (4), pistons (6), piston mount (7) and springs (5) (**these components only**) with an appropriate cleaning agent (e.g., water, alcohol. Do not use organic solvents).
- Let these parts dry completely. Residual moisture in the nose cone assembly may result in a loss of accuracy.

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## Assembling the instrument

- Put the springs on the pistons.
- Observe the colour marking on one side of the piston mount (7) and nose cone assembly (4).

**Important:** These colour markings must be lined up during reassembly.

- Slide the nose cone assembly over the pistons.
- Tighten both screws to secure the piston assembly to the two metal support rods.
- Slide on the manifold housing (3). Turn the closures into place by rotation through 90°.
- Push the handle and the manifold together and turn until they click into place.  
**Important:** Do not push on the ejector cap during this procedure.

## After cleaning:

- Check if pipetting and ejector functions are correct.  
If this is not the case, dismantle the instrument and carefully reassemble it. If necessary, lubricate the pistons, ➤ see "Maintenance", page 9.
- Check for proper fit of pipette tips (page 9).
- Check for correct volume delivery (page 12).

## Autoclaving

**Caution:**  
Only the manifold may be steam sterilized. Eject the tips before autoclaving!

- Remove the manifold from the handle.
- Sterilize the complete manifold in an autoclave. Withstands steam sterilization at 121 °C, 2 bar absolute, for 20 minutes exposure time, according to DIN. It is the user's responsibility to ensure effective sterilization.

## After autoclaving:

- Allow the manifold to cool down completely.
- Push the manifold and the handle together and turn until they click into place.  
**Important:** Do not push on the ejector cap during this procedure.
- Check if pipetting and ejector functions are correct.
- Check for proper fit of pipette tips (page 9).
- Check for correct volume delivery (page 12).

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## Trouble Shooting

Problem	Possible Cause	Corrective Action
Pipette tip leaks or volume too low	Tip does not sit properly  Damaged V-ring (O-ring) on nose cone Piston damaged  Piston seal damaged  Instrument is dirty	Press tip on firmly. Use only suitable high-quality tips. ➤ page 7  Replace V-ring (O-ring) on nose cone. ➤ "Repairs", page 20 Replace piston. ➤ "Repairs", page 20  Replace nose cone assembly. ➤ "Repairs", page 20 Clean instrument. ➤ "Cleaning", page 15
Poor sample aspiration	Nose cone clogged	Clean nose cone. ➤ "Cleaning", page 15
Volume too large	Pipetting key pressed too far down into the blow-out stroke before aspirating reagent	Operate pipette properly. ➤ "Operating the Pipette", page 6

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For available spare parts, see page 27-28.

### Replacing piston springs

- Dismantle the manifold (► "Cleaning", page 15).
- Insert a new spring.
- Reassemble the manifold.

### Replacing pistons

- Dismantle the manifold (► "Cleaning", page 15).
- Roll of the O-ring from the top of the piston.
- Remove the piston from the piston mount.
- Insert a new piston and place the O-ring on the top.
- Reassemble the manifold (► page 16).

### Replacing the nose cone assembly

- Follow procedure under "Cleaning", page 15.

## Replacement of nose cone V-rings (O-rings)

Required materials

- 1 V-Ring removal tool (square)
- 1 V-Ring mounting tool (cone shaped)\*
- Replacement V-rings (O-rings)

\* The mounting tool is only necessary for replacing the V-rings on instruments of the volume ranges 10–100 µl, 20–200 µl and 30–300 µl.

**Note:** Each volume range requires a different sized V-ring (o-ring); ► see "Ordering information", page 27.

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### Mounting a new V-ring (O-ring):

4. Slide the cone-shaped mounting tool over the nose cone. Due to shape of the nose cones, the mounting tool is not needed for instruments of volume ranges 0.5-10 µl, 2-20 µl, 2.5-25 µl or 5-50 µl. The V-ring (O-ring) can be positioned directly onto the nose cone.

**Note:** Moisten mounting tool, nose cone and V-ring (O-ring) with dest. H<sub>2</sub>O in order to facilitate mounting.

5. Place the new V-ring onto the mounting tool. The widest side of the V-ring should face towards nose cone manifold. Push the V-ring (O-ring) over the mounting tool into the groove on the nose cone.

6. Remove the mounting tool. Check correct seat of the V-ring (O-ring). If twisted, turn it several times around the nose cone to restore the V-ring (O-ring) to its correct shape.

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**Important:**  
Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

- Therefore, contact BrandTech Scientific, Inc. and obtain return authorization **before** sending your instrument for service.

**Note:**  
No warranty claims can be considered if repairs other than those described have been made or attempted.

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\* Final test values relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C) of instrument, ambience and dist H<sub>2</sub>O, and with smooth, jerk-free operation. The absolute error margins (µL) are not exceeded at any partial volume.

\*\* Corresponding to the color of the tips.

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Following tolerances are typically (manufacturer's experience) obtained by new instruments under optimized testing conditions (qualified operators and standardized ambience conditions):

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