# Transferpette®

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## Contents

### Seite

		Page
4	Safety Instructions	4
4	Purpose	4
5	Limitations	5
5	Operating Elements	5
6	Handling	6
6	- Use proper pipette tips	6
6	- Volume Setting	6
6	- Hand posture while pipetting	6
7	- Ejecting the tip	7
8	Pipetting	8
8	- Aspirating the sample	8
9	- Dispensing the sample	9
10	Maintenance	10
10	- Lubricating the plunger	10
11	Calibration	11
11	- Fixed-volume models	11
12	- Variable-volume models	12
14	Checking the volume	14
15	Calculation of accuracy/coefficient of variation	15
16	Cleaning	16
19	Autoclaving	19
21	Troubleshooting	22
23	Specifications / Spares	23
26	Accessories	26
27	Return to repair	27
28	Warranty	28

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Read carefully before use!

## Safety Instructions

1. Observe general hazard prevention instructions and safety regulations; e.g. wear protective clothing, goggles and gloves.
2. Observe the reagent manufacturer's information.
3. Read this Operating Manual carefully.
4. Always use the instrument in such a way that neither the user nor any other persons are endangered.
5. Only use manufacturer's original parts.
6. Operate only with tip attached.

## Purpose

Air-displacement pipette for pipetting aqueous solutions of medium density and viscosity.

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
## Operating Limitations

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:

- for liquids incompatible with polypropylene,
- for liquids of a very high vapour pressure,
- in Reverse Mode (ISO 8655/2).

Viscous and wetting liquids may impair the accuracy of the volume. The same applies to liquids of a temperature deviating from the ambient temperature by more than  $\pm 5^{\circ}\text{C}$  (permissible operating range: 15 - 40 °C).

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

## Handling

- Use proper pipette tips  
Pipette performance can only be guaranteed if high-quality pipette tips are used. We consequently recommend the use of PLASTIBRAND® pipette tips.

### Fitting the tip

- turn them lightly in the direction of the arrow

#### Note:

If the tip is turned in the other direction, the screw connection of the shaft tip may be loosened, with air leaks resulting. The shaft tip must be tightened at the latest when a red ring becomes visible on it.

#### Ⓐ Volume setting

#### Ⓑ Hand posture while pipetting


- Put your thumb across the pipetting key
  - unlike with conventional pipettes.

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

## Pipetting

### Take up sample

- Rinse the tip once with the sample liquid.
- 1. Press the lateral pipetting key to **the first stop (A)**.
- 2. Immerse the pipette tip 2 - 3 millimeters into the sample.
- 3. Let the pipetting key slide back **slowly**.

#### Note:

Leave the tip immersed in the liquid for about another second to avoid the intake of air.

- 4. Lightly wipe the tip against the wall of the vessel.

#### Note:

Don't lay the instrument horizontal when the tip is full, or liquid may penetrate inside the instrument.

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
## Discharging the sample


Place the pipette tip against the wall of the vessel.

- Press the pipetting key down to the stop (A) and hold it down.

### Note:

For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- The blow-out stroke empties the tip completely: Press down to the stop (B).

### Note:


When working with serum and liquids of high viscosity, low surface tension or at low volumes: Rinse the tip with the reagent contained in the vessel to improve the accuracy.

- Wipe the pipette tip against the wall of the vessel.
- Let the pipetting key slide back.

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## Maintenance

### In regular intervals:

1. Check firm seat of pipette shaft cone, tighten if necessary.
2. Examine piston and seal for contamination (see page 16 “Cleaning”).
3. Check volume.

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## Lubricating the plunger

Very thinly lubricate the plunger:

- after cleaning,
- if motion is uneven,
- after autoclaving.

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
## Adjustment

The instrument is permanently adjusted for aqueous solutions. It can be adjusted for solutions of other density and viscosity.

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
1/250 vom Nennvolumen

#### Note:

Maximum permissible adjustment is 3 full turns.

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

## Variable-volume models


1. Check the volume, determine actual value (see page 14, "Checking the volume").
2. Release lateral closures by turning into position ① .
3. Pull off ejector cap (A).

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4. Use a paper clip to extract the red distance plate (B) from the top.
5. Pull out volume setting knob (C) to the stop in the direction of its axis (approx. 2 mm).

6. Set to the previously measured actual value. Push back volume setting knob to the stop (approx. 2 mm), without turning it.
7. Replace the red distance plate.
8. Mount ejector cap (A).
9. Secure ejector cap by turning the lateral closures into position  $\ominus$ .

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## Checking the volume

The gravimetric volume test of the pipette requires the following steps:  
(Testing procedures are described e.g., in ISO 8655 Part 6, and DIN 12650 Part 6).

Adjust volume to nominal capacity.

### Conditioning the pipette:

To condition the pipette before testing, mount a new pipette tip and pipette five times with testing liquid ( $H_2O$  dist.). Discard the pipette tip.

### Carrying out the test (observe instructions by the balance manufacturer):

- Mount a new pipette tip and pre-rinse once with testing liquid.
- Take up testing liquid and pipette into the weighing vessel.
- Weigh the pipetted quantity with an analytical balance.
- Calculate the volume, taking the temperature into account.
- Carry out at least 10 pipetting and weighing operations.

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### Variationskoeffizient

$$VK\% = \frac{100 s}{\bar{V}}$$

## Calculation of accuracy (A%) and coefficient of variation (CV%):

Accuracy (A%) and coefficient of variation (CV%) are calculated according to the equations used in statistical quality control.

### Calculations (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n} \quad x_i = \text{results of weighings} \\ n = \text{number of weighings}$$

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
# Cleaning

## Note:

After cleaning, always check volume (see page 14).

## Pipettes up to 200 µl

### (Grey and yellow colour code)

1. Eject the tip.
2. Unscrew pipette shaft (C) from the handle.
3. Remove plunger (A) and plunger spring (B) from the handle.
4. Unscrew pipette shaft cone (D) from the pipette shaft (C).
5. Clean parts shown on drawing with soap solution or isopropyl alcohol. Afterwards rinse with distilled water.
6. Dry the parts (max. 120 °C).
7. Lubricate plunger very thinly (“Maintenance”, page 10).
8. Reassemble parts after they have cooled down.

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
## Pipettes of 200 µl upwards

1. Eject the tip.
2. Unscrew pipette shaft (C) from handle.
3. Unscrew retention bushing (F) from pipette shaft (C).
4. Remove plunger (A) and plunger spring with seal (B), (E).

### Note:

Do not disassemble parts (B) and (E) any further. Take care not to deform the seal; clean it only with a soft lint-free cloth.

5. Unscrew pipette shaft cone (D) from the pipette shaft (C).
6. Clean parts shown on drawing with soap solution or isopropyl alcohol. Afterwards rinse with distilled water.
7. Dry the parts (max. 120 °C).
8. Lubricate plunger ("Maintenance", page 10).
9. Reassemble parts after they have cooled down.

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– continued –

### Note:

Do not overtighten the retention bushing (F). Tighten only by hand.

10. After assembling, operate the pipetting lever five times to centre the seal.

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
## Autoclaving

The pipette parts shown here withstand steam sterilization at 121 °C (2 bar absolute) for 20 min. exposure time ( $t_e$ ), according to DIN.

### Attention:

It is user's responsibility to ensure effective autoclaving.

### After each autoclaving:

- Check volume (see page 14).

### Note:

Only autoclave parts that have been cleaned.

## Pipettes up to 200 µl

### (Grey and yellow colour code)

1. Eject the tip.
2. Cautiously unscrew pipette shaft (C) from handle.
3. Remove plunger (A) and spring (B) from the handle.

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
## Pipettes up to 200 µl

- continued -

4. Sterilize the complete shaft, plunger and plunger spring in the autoclave.
5. Allow the parts to cool completely and assemble in the reverse order.

## Pipettes of 200 µl upwards

1. Eject the tip.
2. Cautiously unscrew pipette shaft (C) from handle.
3. Autoclave the complete shaft (C) in the autoclave without having it disassembled.
4. Allow to dry and to cool completely.
5. Screw pipette shaft into handle.

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## Troubleshooting

Trouble	Possible Causes	Action to be taken
<ul style="list-style-type: none"><li>- Tip dripping (instrument leaks)</li><li>- Aspirated volume too low</li></ul>	<ul style="list-style-type: none"><li>- Unsuitable tip</li><li>- Tip not seated tightly</li><li>- Pipette shaft cone loose</li><li>- Seal contaminated</li><li>- Seal damaged</li><li>- Non-aqueous solutions pipetted</li><li>- Piston damaged</li></ul>	<ul style="list-style-type: none"><li>- Only use high-quality tips</li><li>- Press tip on firmly</li><li>- Tighten clockwise see "Handling", page 6</li><li>see "Cleaning", page 16</li><li>see "Spares", page 25</li><li>see "Adjustment", page 11</li><li>see "Spares", page 25</li></ul>
<ul style="list-style-type: none"><li>- Slow suction or none</li></ul>	<ul style="list-style-type: none"><li>- Shaft clogged</li></ul>	<ul style="list-style-type: none"><li>- Clean shaft see "Cleaning", page 16</li></ul>
<ul style="list-style-type: none"><li>- Volume too high</li></ul>	<ul style="list-style-type: none"><li>- Pipetting key pressed too far (into the blow-out position) before sample take up.</li></ul>	<ul style="list-style-type: none"><li>- Operate properly</li><li>see "Pipetting", page 8</li></ul>

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# Specifications/Spares

(Spare parts can be fitted by the user)

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\* Final test values relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C) of instrument, ambience and dist H<sub>2</sub>O, and with smooth, jerk-free operation. The absolute error margins (µl) are not exceeded at any partial volume.


\*\* Corresponding to the color of the tips.

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Following tolerances are typically (manufacturer's experience) obtained by new instruments under optimized testing conditions (qualified operators and standardized ambience conditions):

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\*\*\* Pipette shaft complete with parts (A) to (F)

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## Accessories

<b>Rack</b> 3 instruments,	Cat. No. 7032 03
<b>Rack</b> 2 x 3 instruments,	Cat. No. 7032 08
<b>Wall support</b> for 3 instr.,	Cat. No. 7032 10
<b>Silicon oil</b>	
V 140 000 mm <sup>2</sup> /s,	Cat. No. 7032 05
<b>Spring</b> (C)	Cat. No. 9884
<b>Retention bushing</b> (A)	Cat. No. 6530 92

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**Note:**

ISO 9000 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend to check the volume every 3 - 12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

The detailed testing instruction can be downloaded on [www.brand.de](http://www.brand.de) or can be ordered under Cat. No. 9976 19 (english) in printed form.

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## Return for Repair

### Attention:

Transporting of hazardous materials without a permit is a violation of federal law.

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BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

- Therefore, contact BrandTech Scientific, Inc. and obtain return authorisation **before** sending your instrument for service.

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## **Warranty, 1 year**

As provided by law, any and all warranties are null and void if the product has been misused, modified or repaired by unauthorized personnel, if the defects are caused by negligence (instruction manual, maintenance) or by normal wear and tear. Use only original manufacturer's accessory/spare parts.

For conditions and extent of warranty refer to our General Conditions of Sale ([www.brandtech.com](http://www.brandtech.com)).

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### **© Spitze abwerfen**
BrandTech Scientific, Inc. wird keine Geräte annehmen, die nicht ordnungsgemäß gereinigt und dekontaminiert wurden.