Transferette®-8/-12

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Read carefully before you start!

Safety Instructions

- Observe general hazard prevention instructions and safety regulations, e. g., wear protective clothing, goggles and gloves.
- Observe the reagent manufacturer's information.
- Read this Operating Manual carefully.
- Always use the instrument in such a way that neither the user nor any other persons are endangered.
- Only use manufacturer's original parts.
- Operate only with tips attached.

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Purpose

The multichannel pipette is an air displacement system for pipetting aqueous solutions with medium density and viscosity into microtiter plates.

Operating Limitations

The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used for:

- liquids incompatible with polypropylene
- liquids of a very high vapour pressure.

Viscous and wetting liquids may impair the accuracy of the volume. The same applies to liquids of a temperature deviating from the ambient temperature by more than $\pm 5^{\circ}\text{C}$ (permissible operating range: 15 - 40 °C).

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Handling

Use proper pipette tips

Pipette performance can only be guaranteed if high-quality pipette tips are used. We consequently recommend the use of PLASTIBRAND® pipette tips.

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
Adjusting the manifold

The **manifold** can be rotated counterclockwise to the most convenient pipetting position.

Counterclockwise rotation
for ease of handling

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
Note: The manifold must be completely mounted to the handle for accurate pipetting. During the counterclockwise rotation, a clicking noise must be heard. Turning in the opposite direction will loosen the manifold.

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Hand posture while pipetting

Securely attach tips to nose cone assembly. Put your thumb **across** the pipetting key – unlike with conventional pipettes.

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

Pipetting

Aspirating the sample

- Rinse the tips once with the sample liquid (at normal pipetting speed).
- Press the pipetting key on the side to the **first stop (A)**.
- Dip the pipette tips 2-3 millimeters into the sample.

Note: Leave the tips in the liquid for about one second longer to avoid the intake of air.

- Let the pipetting key slide back **slowly**.
- Lightly wipe the tips against the wall of the vessel.


Note: Don't lay the instrument horizontally when the tips are full, or liquid may flow inside the instrument.

Dispensing the sample

Lightly place the pipette tips against the wall of the vessel.

- Press the pipetting key down to the **stop (A)** and hold it down.

Note: For serums and liquids of high viscosity or low surface tension, wait about 3 seconds. This helps to improve accuracy.

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After use:

Always store the instrument upright on the **stand** supplied.

Note: To reduce strain on the sealing rings (► "V-rings", page 22), always store the instrument without the tips attached.

Ejecting the tips (A)


Press the ejector cap down firmly (see arrow).
Tips that are jammed on too tightly may need to be pulled off separately by hand.

Volume setting (B)

Set the desired volume by turning the volume setting knob (► "Volume setting knob", page 4).

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

Maintenance

Check for tightness

After drawing in the sample:
Hold down the instrument vertically for about 10 seconds. If a drop forms at one of the tips,

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
Lubrication

After autoclaving the manifold several times, it is necessary to lubricate the pistons to maintain smooth operation:

Unscrew nose cones, push pistons downward, and lubricate (► "Repairs", page 18, 19).

To lubricate the pistons, use only the silicone grease supplied with the instrument.

The silicone grease is also available as accessory (► "Ordering Data", page 31).

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
Adjustment

The instrument is permanently calibrated for aqueous solutions. It can be adjusted for solutions of other density and viscosity.

1. Check the volume, determine actual value (► page 13).
2. Release lateral closures by turning into position ①.
3. Remove ejector cap (A).
4. Use a paper clip to extract the red security plate (B) from the top.
5. Pull out volume setting knob (C) to the stop in the direction of its axis (approx. 2 mm).

6. Set to the previously measured actual value. Push back volume setting knob to the stop (approx. 2 mm), without turning it.
7. Replace the red security plate (across the axis).
8. Mount ejector cap (A).
9. Secure ejector cap by turning the lateral closures into position ②.

Checking the volume

The instrument should be calibrated at least once a year. A more frequent interval may be dictated by individual lab rules, applications or conditions. The gravimetric volume test of the pipette requires the following steps (page 12):
(Testing procedures are described e.g., in ISO 8655 Part 6 ► "Testing Instruction", page 34.)

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– Adjust volume to nominal capacity.

Conditioning the pipette:

– To condition the pipette before testing, mount new pipette tips and pipette five times with testing liquid (H_2O dist.). Discard the pipette tips.

Carrying out the test (observe instructions by the balance manufacturer):

- Mount new pipette tips and pre-rinse once with testing liquid.
- Take up testing liquid and pipette into the weighing vessel.
- Weigh the pipetted quantity with an analytical balance.
- Calculate the volume, taking the temperature into account.
- The number of tests per channel should be selected according to the required accuracy. We recommend 5-10 tests per channel in 3 volume ranges (100 %, 50 %, 10 %).

Calculation of accuracy (A%) and coefficient of variation (CV%): A and CV are calculated according to the equations used in statistical quality control.

Calculations (for the nominal volume)

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

Components

- ① Handle
- ② Manifold
- ③ Manifold housing
- ④ Nose cone assembly
- ⑤ Individual nose cone
- ⑥ Piston springs
- ⑦ Piston
- ⑧ Piston support bar
- ⑨ Manifold housing cover
- ⑩ Seal
- ⑪ Mounting tool
- ⑫ Central guide rod

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
Cleaning

Dismantling the instrument


- Be careful:** Work over a tray to collect any part that may fall out during dismantling.
- Eject the tips.
 - Set the instrument to nominal volume.
 - Remove the manifold (②) from the handle (①).
 - Turn both closures of the manifold housing cover (⑨) with a coin (rotation through 90°) and slide off the housing.
 - **Remove both outer screws (Phillips). Do not loosen the central guide rod (⑫)!**
 - Remove the nose cone assembly (④), piston springs (⑥), and piston support bar (⑧).
 - Unscrew single nose cones with seals from the nose cone support bar (► "Repairs", page 18).

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Cleaning the components

- Clean single nose cones (⑤), pistons (⑦), piston support bar (⑧) and piston springs (⑥) (**these components only**) with an appropriate cleaning agent (e.g., water, alcohol. Do not use organic solvents!).
- Let these parts dry completely. Residual moisture in the nose cone assembly may result in a loss of accuracy.
- Lubricate the pistons lightly (only with the supplied silicone grease).

Important: To lubricate the central guide rod, use only the recommended fluorostatic grease (► page 31)!

Assembling the instrument

- Reassemble the nose cone assembly without cones (④), springs (⑥), piston support bar (⑧) and manifold housing cover (⑨). Tighten all screws.
- Set volume to minimum.
- Push the seals on the pistons with the flat side facing upward (direction of the handle) (► Fig. 4, page 19); use the mounting tool to mount and tighten the single nose cones clockwise.
- Slide on the manifold housing (③). Turn the closures into place by rotation through 90°.
- Push the handle and the manifold together and turn until they click into place.

Important: Do not press on the ejector cap during this procedure.

After cleaning:

- Check if pipetting and ejector functions are correct. If this is not the case, dismantle the instrument and carefully reassemble it. If necessary, lubricate the pistons (► "Maintenance", page 9).
- Check for proper fit of tips.
- Check for correct volume delivery (► page 11).

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



Repairs

- "Accessories and spare parts", pages 31-33.

Replacing of single nose cones and seals

Removing a single nose cone

- Set the instrument to nominal volume.
- Turn both closures of the manifold housing cover (⑨) with a coin (rotation through 90°) and slide off the housing.
- Push the mounting tool (⑩) on a single nose cone (Fig. 1).
- Turn counterclockwise to loosen and remove the nose cone (Fig. 2).

Note: The seal will either remain inside the nose cone, or will stay on the piston after the nose cone is removed. By pressing the pipetting key it can easily be removed and cleaned (Fig. 3).

- Set volume to minimum.
- Push down the pipetting key completely.
- If required, lubricate the piston lightly with the **supplied** silicone grease.
- Push the seal on the piston with its flat side facing upward (direction of the handle) (Fig. 4).
- Use the mounting tool to mount and tighten the new nose cone clockwise.
- Reassemble the manifold (► page 16).
- Check instrument for tightness (► "Check for tightness", page 9).

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



Replacing the pistons

- Unscrew the manifold from the handle.
- Turn both closures of the manifold housing cover (⑨) with a coin (rotation through 90°) and slide off the housing.
- **Remove both outer screws (Phillips)** which hold the nose cone assembly (④) to the metal support rods of the manifold housing cover (⑨).
Do not loosen the central guide rod (► page 14)!

- Push the mounting tool (11) on a piston (7).
- Turn counterclockwise to loosen and remove the piston (7) along with the spring (13).
- Insert a new piston.
- Use the mounting tool to mount and tighten the new piston.
- Reassemble the instrument (► page 16) and check for tightness (► page 9).

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Mounting a new V-ring (O-ring):

4. Only for instruments of size 2.5-25 µl, 5-50 µl, 10-100 µl, 20-200 µl and 30-300 µl: Slide the cone-shaped mounting tool over the nose cone.

For instruments of size 0.5-10 µl and 2-20 µl, no mounting tool is required; the O-ring can be placed directly onto the nose cone.

Note: Moisten mounting tool, nose cone and V-ring (O-ring) with dist. H₂O in order to facilitate mounting.

5. Place the new V-ring onto the mounting tool. The widest side of the V-ring should face towards nose cone manifold. Push the V-ring (O-ring) over the mounting tool into the groove on the nose cone.

6. Remove the mounting tool. Check correct seat of the V-ring (O-ring). If twisted, turn it several times around the nose cone to restore the V-ring (O-ring) to its correct shape.

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

After any repair:

- Check if pipetting and ejector functions are correct.
- Check for proper fit of the tips.
- Check for correct volume delivery (► page 11).

After every change of the complete manifold:

- Check for correct volume delivery (► page 11). Readjust instrument if necessary.

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Autoclaving

Caution: Only the manifold may be steam sterilized. Eject the tips before autoclaving!

- Remove the manifold from the handle.
- Sterilize the complete manifold in an autoclave. Withstands steam sterilization at 121 °C, 2 bar absolute, for 20 minutes exposure time, according to DIN. It is the user's responsibility to ensure effective sterilization and correct autoclave temperature.

After autoclaving:

- Allow the manifold to cool down completely.
- Push the manifold and the handle together and turn until they click into place. **Important:** Do not press the ejector cap during this procedure.
- If necessary, tighten the nose cones with a sterile pipette tip of proper size.
- Check for smooth function of pipetting and ejector keys.
- Check for proper fit of pipette tips.
- Check for correct volume delivery (► page 11).

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Trouble Shooting

Problem	Possible Cause	Corrective Action
Pipette tip leaks or volume too low.	Tip does not sit properly.	Press tip on firmly. Use only suitable high-quality tips, page 4.
	Damaged V-ring (O-ring) on single nose cone.	Replace V-ring (O-ring) on nose cone, "Repairs", pages 22-24.
	Piston damaged.	Replace piston "Repairs", pages 20, 21.
	Seal dirty, damaged.	"Repairs", pages 18, 19.
	Instrument is dirty.	Clean instrument "Cleaning", page 15.
Poor sample aspiration.	Single nose cone clogged.	Clean nose cone. Remove cones only; do not disassemble instrument! "Cleaning", page 15.
Volume too large.	Pipetting key pressed too far down into the blow-out stroke before aspirating reagent.	Operate pipette properly "Operating the Pipette", page 6.
Pipetting key moving stiffly.	Key movement reduced.	Lubricate piston and central guide rod with the appropriate grease (pages 9, 31).

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* These final test values refer to the delivered nominal volume, instrument and liquid (H_2O dist.) at equilibrium with ambient temperature (20 °C), and smooth, jerk-free operation. The absolute error (given in µl) is not exceeded at any partial volume.

** Corresponding to the colour of the tips.

When pipetting liquids with different viscosity and/or flow properties than H_2O dist., deviations from the values stated (DIN 12650-6) can occur.

When calibrating according to the rules for monitoring of measuring instruments, we recommend the use of "general practice tolerances". However, these should not be higher than double the tolerances as stated in DIN 12650-6.

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Note:

ISO 9000 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

The detailed testing instruction can be downloaded from www.brand.de or can be ordered under Cat. No. 9976 23 (English) in printed form.

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Return for Repair

Important: Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

► Therefore, contact BrandTech Scientific, Inc. and obtain return authorization **before** sending your instrument for service.

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Warranty, 1 year

As provided by law, any and all warranties are null and void if the product has been misused, modified or repaired by unauthorized personnel, if the defects are caused by negligence (instruction manual, maintenance) or by normal wear and tear. Use only original manufacturer's accessory/spare parts.

For conditions and extent of warranty refer to our General Conditions of Sale (www.brandtech.com).

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BrandTech Scientific, Inc. wird keine Geräte annehmen, die nicht ordnungsgemäß gereinigt und dekontaminiert wurden.
Brandtech Scientific, Inc. n'accepte que des appareils soigneusement nettoyés et décontaminés.
BrandTech Scientific, Inc. no acepta apartos que no hayan sido debidamente limpiados y descontaminados.