# Macro Transferpettor™

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#### Please read carefully before use!

# Safety Instructions

This instrument may be involved with hazardous materials, operations, and equipment. This manual does not purport to address all of the safety problems associated with its use. It is the responsibility of whoever uses this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

- 1. Observe the general instructions for hazard prevention and safety regulations. E.g. wear protective clothing, goggles and protective gloves.
- 2. Observe the reagent manufacturers' information.
- 3. Read and follow this Operating Manual carefully step by step.
- 4. Only use the instrument for pipetting and only within the resistance limits of the materials. Observe the operating limitations.
- Always use the instrument in such a way that neither the user nor any other person is endangered.

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- 6. Dangerous furning chemicals must be pipetted in a fume hood.
- 7. Never use force.
- 8. Only use original accessories and original spare parts.

# **Operating Limitations**


Liquids which attack PP or PE are not to be pipetted.

The instrument does not withstand steam sterilization.

Crystallization of concentrated saline solutions and other crystallizing liquids may damage the seal.

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## **Purpose**

The instrument is designed for pipetting liquids. It employs the principle of positive displacement. The piston (2) (= Transferpettor-seal) draws the liquid directly, without air interface, into the tip (1) (= Transferpettor-cap).

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## Range for use

Suitable for media with:

- Density up to 13.6 g/cm<sup>3</sup>
- Viscosity up to 50 000 mm<sup>2</sup>/s
- Vapour pressure up to 500 mbar

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## Setting the volume

- Set the lever (1) to  $\widehat{\mathbf{\Omega}}$  .
- Turn the pipetting key (2) until the required volume is displayed (3).
- Set the lever (1) to  $\widehat{\Omega}$ .

Recommended working range: 20 to 100 % of the nominal volume. At volume settings below 20 % precision can no longer be guaranteed.

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
## **Pipetting**


The instrument is immediately ready for use, as supplied.

## To draw in liquid

- Push the pipetting key fully down.

- Immerse the tip (= Transferpettor-cap) into the liquid (2 to 3 mm).
- Slowly release the pipetting key.

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
 Carefully wipe the outside of the tip with a lint-free cloth or swab. Take care not to touch the tip orifice because this would cause liquid to be drawn from the tip.

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
- To dispense the liquid
- Place the tip (= Transferpettor-cap) against the wall of the vessel.
- Push the pipetting key down.

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
- Lightly wipe the tip against the wall of the vessel.
- Release the pipetting key.

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## Daily:

# **Checking for Tightness**

- Set the nominal volume.
- Aspirate in liquid.
- Hold the instrument upright for 10 sec, (tip down).
- If a drop forms, > "Trouble Shooting", page 22.

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
## Maintenance

Maintenance is not required, if the instrument is properly used.

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Changing the Tip
(= Transferpettor-cap)

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

 Turn the tip until the arrow on the instrument and the arrow on the tip are in alignment.

Do not touch the contaminated area of the

tip (= Transferpettor-cap) and the seal.

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

## Attaching the new tip

- Slide the tip (= Transferpettor-cap) over the seal.
- Align the arrow on the instrument and the arrow on the tip.
- Slide the tip into the mounting, keep it slightly pressed into it.
- To position the tip exactly, press the pipetting key fully down and
- fix the tip by turning it to the right or to the left (max. 45°).
- Let the pipetting key slide back.

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#### Note:

If the pipetting key is not fully pushed down when the tip (= Transferpettor-cap) is installed, the tip will be seated too far into the tip mount. As a consequence the instrument will deliver a lower volume than indicated.

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
### Removing the old seal

#### Attention!

Do not touch the contaminated area of the tip (= Transferpettor-cap) and the seal. Wear protective gloves.

- Remove the tip (= Transferpettor-cap).
   "Changing the Tip", page 11.
- Hold the instrument in a vertical position.
- Place the mounting tool on the seal (with the conical opening upwards).
- Push the pipetting key fully down.
- Using the mounting tool, vigorously pull the old seal off the piston rod.
- Let the pipetting key slide back.

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### Attaching the new seal

- Place the new seal into the conical opening of the mounting block.
- Press the piston rod into the seal as far as it will go.
- Check the seal for tight seating.
- Slide the tip (= Transferpettor-cap) on.
  - > "Attaching the new tip", page 12.

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#### Removing the old piston rod

#### Attention!

Do not touch the contaminated area of the tip (= Transferpettor-cap) and the seal. Wear protective gloves.

- Remove the tip (= Transferpettor-cap)."Changing the Tip", page 11.
- Push the pipetting key fully down.
- Turn the piston rod locking device anticlockwise to unscrew the piston rod from the instrument.
- Disconnect the piston rod locking device from the piston rod.

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## Attaching the new piston rod

- Screw the piston rod locking device on the piston rod to lock position.
- Then screw the piston rod into the instrument to lock position.
- Push the pipetting key down and finger-tighten the piston rod locking device.
- Attach a new seal, > page 14.
- Attach a new tip, > page 12.

# Cleaning

- Dismantle the instrument.
  - > "Changing the Piston Rod", page 15.
- Clean the piston rod.
- Allow to dry completely.
- Reassemble the instrument.
  - > "To fit the new piston rod", page 16.

## Checking the Instrument

after every change of the tip (= Transferpettor-cap), seal, or piston rod

- Check that the pipetting key can be pushed fully down to the end plate (1).
- Ascertain that the tip fits tightly."Changing the Tip", page 11.
- Check the tightness, > page 9.
- Check the volume, > page 18.

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#### Note:

Temperatures exceeding 40 °C and storage periods of over three months may cause a deformation of tip (= Transferpettor-cap) and seal, adversely affecting the tightness of the seal and compromising instrument performance.

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$$V \% = \frac{100 \text{ s}}{\overline{V}}$$

# Checking the Volume

**Note:** The detailed Test Instructiones are available under Cat. No. 9976 27.

The accuracy and coefficient of variation of the instrument are determined gravimetrically as follows:

- Set to the nominal volume.
- Pipette distilled H<sub>2</sub>O.
- Weigh the pipetted quantity on an analytical balance
- Calculate the pipetted volume taking the temperature into account.
- Perform at least 5 pipetting and weighing operations.
- Calculate the accuracy (A%) and the coefficient of variation (CV%) by means of the formulae used in statistical quality control. The proceeding is described e.g. in ISO DIS 8655/2+/3. Observe the operating manual of the balance manufacturer and the corresponding standards.

## Calculations (for the nominal volume)

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# **Trouble Shooting**

Problem	Possible Cause	Corrective Action			
Pipetted volume too low/ Pipetting key cannot be pushed fully down to the end plate.	The tip has been installed too far into the tip mount. Piston rod is loose.	Remove the tip (= Transferpettor-cap). > "Changing the Tip", page 11. Tighten the piston rod locking device. > "To fit the new piston rod", page 16. Attaching the tip, > page 12.			
Pipette leaks.	The tip does not fit properly.	Reinstall the tip (= Transferpett <u>or</u> -cap) and slide it on again. > "Changing the Tip", page 11.			
	Damaged seal.	Change the seal, > page 13.			
Air bubbles in the liquid drawn in.	The liquid was drawn in too quickly.  The tip does not fit properly.  Damaged seal.	Let the pipetting key slide back slowly.  Reinstall the tip (= Transferpettor-cap) and slide it on again.  > "Changing the Tip", page 11.  Change the seal, > page 13.			
Liquid above the seal.	Damaged tip or seal.	Change tip (= Transferpettor-cap) or seal. > "Changing the Tip", page 11. > "Changing the Seal", page 13.			
Liquid remaining in the tip.	Tip not properly seated. Piston rod is loose.	Remove the tip (= Transferpettor-cap).  > "Changing the Tip", page 11.  Tighten the piston rod locking device.  > "To fit the new piston rod", page 16.  Fit the tip, > page 12.			
	Damaged seal.	Change the seal, > page 13.			

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- These values refer to the delivered nominal volume, H<sub>2</sub>O dist., instrument, and medium at equilibrium with ambient temperature (20 °C), and smooth, jerkfree operation.
- 2) Pack of 10.
- 3) Pack of 3.
- 4) Pack of 1.

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Best.-Nr. 7028 90

## Transferpettor Station

To accommodate two instruments of 0.5 to 10 ml, complete with their accessories.

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# Return for Repair

#### Important:

Transporting of hazardous materials without a permit is a violation of federal law.

BrandTech Scientific, Inc. will not accept instruments that are not appropriately cleaned and decontaminated.

➤ Therefore, contact BrandTech Scientific, Inc. and obtain return authorisation <u>before</u> sending your instrument for service.

#### Note:

No warranty claims can be considered if repairs other than those described have been made or attempted.

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BrandTech Scientific, Inc. wird keine Geräte annehmen, die nicht ordnungsgemäß gereinigt und dekontaminiert wurden.