


# epT.I.P.S.<sup>®</sup> Replacement Options



022491903 – 10 µL S Tip (Eppendorf Quality)



022491806 – 10 µL S Tip (PCR Clean)

Keep Purity

Keep Rack



022491504 – 10 µL S Tip Reload

<b>Size</b>	Identical
<b>Format/Fit</b>	Reload
<b>Purity</b>	Identical



0030071557 – 10 µL M Tip (Sterile)

<b>Size</b>	Identical*
<b>Format/Fit</b>	Identical
<b>Purity</b>	Sterile (Upgrade)

\* Except 10 S and 0.5–20 L Tips

Keep Purity

Keep Rack



022491709 – 10 µL S Tip Reload

<b>Size</b>	Identical
<b>Format/Fit</b>	Reload
<b>Purity</b>	Identical



0030078500 – 10 µL S Tip (Filter)

<b>Size</b>	Identical
<b>Format/Fit</b>	Identical
<b>Purity</b>	PCR-filter (Upgrade)

# epT.I.P.S.<sup>®</sup> Replacement Options

## Ordering information

epT.I.P.S. Reloads Eppendorf Quality	Old Order No. for Racks	New Order No. (Same Purity Reload)	New Order No. (Same Rack Purity Upgrade)
0.1–10 µL, 34 mm, dark gray, 960 tips (10 trays x 96 tips)	022491903	022491504	0030071557
0.1–20 µL, 40 mm, medium gray, 960 tips (10 trays x 96 tips)	022491911	022491512	0030071557
0.5–20 µL L, 46 mm, light gray, 960 tips (10 trays x 96 tips)	022491920	022491521	
2–200 µL, 53 mm, yellow, 960 tips (10 trays x 96 tips)	022491938	022491539	0030075234
20–300 µL, 55 mm, orange, 960 tips (10 trays x 96 tips)	022491946	022491547	0030071573
50–1,000 µL, 71 mm, blue, 960 tips (10 trays x 96 tips)	022491954	022491555	0030071581
50–1,250 µL, 76 mm, green, 960 tips (10 trays x 96 tips)	022491962	022491563	0030071590
0.25–2.5 mL, 115 mm, red, 480 tips (10 trays x 48 tips)	022491971	022491571	0030075285
epT.I.P.S. Reloads PCR clean	Old Order No. for Racks	New Order No. (Same Purity Reload)	New Order No. (Same Rack Purity Upgrade)
0.1–10 µL, 34 mm, dark gray, 960 tips (10 trays x 96 tips)	022491806	022491709	0030078500
0.1–20 µL, 40 mm, medium gray, 960 tips (10 trays x 96 tips)	022491814	022491717	0030078519
0.5–20 µL L, 46 mm, light gray, 960 tips (10 trays x 96 tips)	022491822	022491725	0030078527
2–200 µL, 53 mm, yellow, 960 tips (10 trays x 96 tips)	022491831	022491733	0030078551
20–300 µL, 55 mm, orange, 960 tips (10 trays x 96 tips)	022491849	022491741	0030078560
50–1,000 µL, 71 mm, blue, 960 tips (10 trays x 96 tips)	022491857	022491750	0030078578
50–1,250 µL, 76 mm, green, 960 tips (10 trays x 96 tips)	022491865	022491768	0030078594
0.25–2.5 mL, 115 mm, red, 480 tips (10 trays x 48 tips)	022491873	022491776	0030078586

Eppendorf North America, Inc.  
 Phone: 800-645-3050  
 Email: info@eppendorf.com

Eppendorf Canada Ltd.  
 Phone: 800-263-8715  
 Email: canada@eppendorf.com

[www.eppendorf.com](http://www.eppendorf.com)