

Glass Tip Pipette

Nichipet ECO

- The accuracy and precision of a micropipette with the stability and reliability of glass tips.
- The glass tip enables bias free liquid handling for endocrine disruptor's analysis.
- Four different sizes covering a wide volume range. (20 μ L - 10000 μ L)
- Glass tip is secured to the unit with a double O-ring design.
- Comes with a PTFE nozzle and ceramic piston.



Specifications

Cat. No	Volume range(μ L)	Increment (μ L)	Accuracy(%)	Precision(%)	Usable tips
00-NPC-200	20 - 200	0.2	$\pm 6.0 - \pm 3.0$	$\leq 6.0 - \leq 3.0$	Glass Tip S
00-NPC-1000	100 - 1000	1.0	$\pm 5.0 - \pm 3.0$	$\leq 5.0 - \leq 3.0$	Glass Tip L
00-NPC-5000	1000 - 5000	10.0	$\pm 5.0 - \pm 3.0$	$\leq 5.0 - \leq 3.0$	Glass Tip X
00-NPC-10ML	1000 - 10000	50.0	$\pm 6.0 - \pm 3.0$	$\leq 6.0 - \leq 3.0$	Glass Tip Z

● Glass tips

Cat. No	Volume range(μ L)	Q'ty(box)
00-GLT-S	20 - 200	10
00-GLT-L	100 - 1000	10
00-GLT-X	1000 - 5000	10
00-GLT-Z	1000 - 10000	10

Further research into "Environmental Hormone" is strongly requested to clarify the precise cause and effect on the human system and wild living things. Glass tips rather than standard plastic pipette tips are more suitable for this research.

