Portable Pipet-Aid® XP2

The Evolution of Xcellence Continues

- New Ergonomic Design—The most comfortable pipettor you ever laid a hand on
- New Power Source for Uninterrupted Extended Operation— Can be charged while in use
- New Ultra Quiet Precision Pump—Great control for aspiration or dispensing

The new Drummond Pipet-Aid XP2 advances the industry standards for pipetting **comfort**, **convenience**, and **control**. Building on the years of experience of providing state-of-the-art tools that make your job in the laboratory safer, easier, and more convenient, the Pipet-Aid XP2 delivers more. The ergonomic design is more comfortable to handle. The new power source enables more convenient extended operation and the new more powerful pump provides great control and quieter operations. The unit is supplied complete with a power supply/charger, 4 extra filters, and a removable stand.

Specifications

Size
Weight
Capacity Accepts 1 to 100 mL pipets
Battery Ni-MH (3.6V/700 mAh)

Cat. No.	Description
4-000-501	Pipet-Aid XP2, 110V
4-000-501-A	Pipet-Aid XP2, 240V, Australia
4-000-501-I	Pipet-Aid XP2, 240V, CE, International
4-000-501-UK	Pipet-Aid XP2, 230V, CE, UK

Replacement Parts

rtopiaco.	inonic i di		
Cat. No.		Description	
4-000-520	(1)	Nosepiece Cover	
4-000-522	(2)	Filter Lock	
4-000-523	(3)	Handle Shell	
4-000-524	(4)	Filter Connector	
4-000-051	(5)	Self-locking TC Filter 5/pk	
4-000-051-S		Sterile Self-locking TC Filter 5/pk	
4-000-052		Self-locking TC Filter 50/pk	
4-000-052-S		Sterile Self-locking TC Filter 50/pk	
4-000-002	(6)	Rubber Insert 2/pk	(2)
4-000-525	(7)	Button Assembly (XP2 and XL)	Ĭ
4-000-528	(8)	XP2 and XL Valve Body Complete w/Buttons	$\stackrel{\downarrow}{=}$
4-000-072	(9)	Switch Operator, XP/XL/XP2	3
	(10)	Valve Body Slide (yellow) comes with Valve Body	(3)
4-000-561	(11)	XP2 Pump, Board, Wire Harness, Micro Switch with Connecting Tubing	4
4-000-526	(12)	XP2 Removable Stand	(5)-
4-000-527	(13)	Front Cover Hold-down Clip	
4-000-036	(Not Shown)	Battery Pack (XP2 and XL only)	
4-000-065*	(14)	XP2 and XL Power Supply 100 - 240V, 50/60 Hz *Now replaces: 4-000-066 4-000-066-A 4-000-066-E	6
4-000-530	(15)	XP2 Tissue Culture Nosepiece Kit	











